

















**Analytical Chemistry-1-Lab. for
First Stage of Chemistry
semester 2^{ed}**

**Laboratory Manual of
Semimicro Qualitative Inorganic
Analysis**

**Dr. Wafaa Waleed AL-Qaysi
&
Yasmeen H. Muhamad**

2020-2021






ارشادات مختبرية عامة

- لاخذ مادة صلبة من قنينة استعمل ملعقة زجاجية او بلاستيكية لا تسكب المادة الصلبة فوق ورقة الترشيح او انبوبة اختبار او ورق او كاس زجاجي بينما السوائل فقط هي التي تسكب في الاوعية الزجاجية المناسبة لها
- لا تغمر قطارة او ماصة في قنينة سائل او محلول لتستخدمه في سائل او محلول اخر
- لا تضع ابداء سداد قنينة او محلول فوق الرف او المنضدة تحاشيا من التلوث
- ضع جميع قناني المواد الكيميائية في الاماكن المخصصة لها بعد اخذك ما تحتاج منها دون حملها الى مكان عملك فالطلبة الاخرى هم بحاجة لها ايضا
- لا تعيد الكمية الزائدة من المواد الكيميائية(الصلبة كانت ام سائلة الى قنيتها الاصلية) اعطاها لزميل بحاجتها او اخبر احد مسؤولي المختبر لمعرفة ما يجب عمله للاستفادة منها
- لا تترك ابداء قنينة مواد كيميائية مفتوحة فهناك كثير من المواد الكيميائية التي يتغير تركيبها عند تعرضها للجو او ضوء الشمس او الضوء الاصطناعي اضافة الى احتمال تلوثها و تاكد من ان القنينة قد سدت باحكام عند اخذ جزء من المادة التي تحتويها
- ضع المواد الصلبة الزائدة التي لا يمكن الاستفادة منها واوراق الترشيح وبقايا عيدان الثقاب وغيرها في وعاء المهملات ولا تضعها في الحوض او في مجرى المياه و تاكد قبل رمي المادة الصلبة انها لا تسبب اشعال النار فيما لو تعرضت للهواء او انها تعطي روائح غير مرغوب فيها عند رميها في الحوض او مجاري المياه
- لا تسكب الماء ابداء في حامض الكبريتيك المركز لغرض تخفيف الحامض الطريقة السليمة لتخفيف حامض الكبريتيك المركز هي باضافة الحامض المركز ببطيء و عناية الى الماء في وعاء زجاجي(انبوب اختبار مناسب او اسطوانة مدرجة او كاس زجاجي او قنينة زجاجية الخ) مع التحريك المستمر فاذا اقتربت درجة حرارة المحلول من مرحلة الغليان عند ذلك يجب تبريده بسكب ماء الحنفية على الجدران الخارجية للوعاء و الاستمرار بالتحريك لحين انخفاض درجة الحرارة ثم يباشر باضافة الحامض كالسابق
- احذر عند تسخين مادة صلبة او سائلة في انبوبة اختبار ان تكون فوهة الانبوبة موجهة نحوك او نحو زميل مجاور لك او باتجاه أدوات المختبر وكذلك نفس الحال عند غلي محلول او سائل
- لا تحاول ابداء تذوق اي مادة كيميائية سائلة او صلبة تتعامل معها في المختبر و لا تشرب الماء في اي وعاء زجاجي تستخدمه في عملك داخل المختبر اذكر دائما ان من حولك العديد من المواد السامة
- كن حذرا عندما تريد شم سائل او غاز يتصاعد اثناء تفاعل ما و لا تشم بخار السائل او الغاز المتصاعد من قنينة او انبوبة اختبار او ورق التفاعل بتقريب الفوهة من الانف
- اجراء عملية التسخين داخل الهود(خزانات الغاز) لان عملية التسخين يتضمن تحرير ابخرة او غازات ذات رائحة كريهة او سامة
- ينبغي على الطالب عدم لبس الحلي الذهبية و الفضية وغيرهما اثناء العمل خشية لتعرضها للتلف

Analytical chemistry

Is the science of the characterization and measurement of chemicals and also involve separating, identifying and determining the relative amounts of the components in a sample of matter. Chemical analysis is divided into two types:



Qualitative analysis

The term **Qualitative analysis** involves the method by which the constituents of materials are identified. Three different methods of operation used in **Qualitative** analysis depending upon the amount of sample used. Macro on **Qualitative** analysis. For example; is applied when the solution being analyzed has a volume ranging from 5 to 100 ml or more. And in the case of solid samples, an amount of 1 g or more are often used. The semi micro method involves the use of smaller quantities of materials ranging, in solid, from 1 mg. to 1 g. and in liquid from 1 drop (about 0.05 ml) to 1 ml.

Positively charged ions are termed cations (i.e. cathode ions) because, in an electrolysis cell, they migrate toward the negative electrode. Positively charged inorganic ions or inorganic cations are, sometimes, called metallic ions, and are generally produced by the ionization of inorganic salts and bases.

In this scheme twenty four cations of common occurrence are involved. By taken advantage of the differences in chemical properties, those cations are divided into five major group. Four of the groups are precipitated and separated from each other through group reagents, whereas the fifth is a soluble one.

Group I (Ag^+ , Pb^{2+} , Hg_2^{2+}) cations: produce insoluble chlorides so they can be precipitated with dilute HCl, while all other cations remain in solution.


Group II (Cu^{2+} , Bi^{3+} , Cd^{2+} , Hg^{2+} , As^{3+} , Sb^{3+} , Sn^{4+}) cations: produce very insoluble sulfides so they can be precipitated by low amounts of sulfide ion; this can be achieved by adding an acidic solution of H_2S .

Group III (Al^{3+} , Cr^{3+} , Fe^{3+} , Zn^{2+} , Ni^{2+} , Co^{2+} , Mn^{2+}) cations: produce slightly soluble sulfides so they can be precipitated by relatively high amounts of sulfide ion; this can be achieved by adding a basic solution of H_2S .

Group IV (Mg^{2+} , Ca^{2+} , Sr^{2+} , Ba^{2+}) cations, as well as all of the above groups, produce insoluble carbonates so they can be precipitated by the addition of carbonate once the ions of the first three groups have been removed.

Group V (Na^+ , K^+ , NH_4^+) cations do not precipitate with any of the above reagents.

Group I Analysis



The silver group is composed of three cations namely silver ion Ag^+ , mercury (I), Hg_2^{2+} , and lead (II), Pb^{2+} . Chemistry of the Precipitation of the Group I Cations The precipitating reagent of Group I is a dilute solution of hydrochloric acid, and the ions precipitated are those of silver, Ag^+ , mercury (I), Hg_2^{2+} , and lead (II), Pb^{2+} .

A. Identification of the Silver ion Ag^+

A list of some of the insoluble Silver compounds is given in table 1 below:


Table 1 Insoluble silver compounds with their colors

Compound	color
$\text{Ag}_2\text{C}_2\text{O}_4$	White
* Ag_2CrO_4	Red
Ag_2CO_3	Yellow
AgOH	Brownish black
AgCl	White
AgCN	White
AgCNS	White
AgBr	Pal yellow
AgI	Yellow
Ag_2S	Grey-black

* Decreasing Concentration of Ag^+ Ion

1. Chloride ion solution (use dilute HCl)


When a dilute hydrochloric acid is added to silver ions solution. a white precipitate of silver chloride, AgCl is formed.



- a) AgCl chemical decomposition when exposed to ultraviolet or to a strong sun light (distinction from PbCl_2 and Hg_2Cl_2). This method or identification is very sensitive with U.V. light. Place a drop of the test solution on a filter paper, followed by a drop of dilute HCl. Expose the filter paper containing the precipitation to U.V. light. Appearance of a purple-gray or grayish-black spot indicates the silver ion.
- b) AgCl is soluble in aqueous ammonia (ammonia water ammonium hydroxide) due to the formation of the silver ammonia complex $\text{Ag}(\text{NH}_3)_2^+$




$\text{Ag}(\text{NH}_3)_2^+$ ion is destroyed upon acidifying the solution with nitric acid reprecipitating AgCl .




2. 8-quinolinol solution (8-hydroxy quinolone or oxime).

When this reagent is added to a solution containing silver ions. a yellow precipitate of silver 8- quinolinolate, $\text{AgC}_9\text{H}_6\text{NO}$, is formed, in the present of ammonium acetate.




3. Iodide ion solution(use KI solution)


Iodide ion precipitates the silver ion as yellow silver iodide, AgI .



AgI is soluble in both solutions of potassium cyanide.




And sodium thiosulfate



4. Chromate ion solution(use K_2CrO_4 solution)

Silver ion and chromate ion combine to form a red precipitate of silver chromate, Ag_2CrO_4 , which is insoluble in dilute acetic acid but soluble in nitric acid and in aqueous ammonia.



B. Identification of the Lead ion Pb^{2+}

A list of some of the insoluble mercuric compounds is given in table 2 below:


Table 2 Insoluble lead compounds with their colors

Compound	color
PbCl_2	White
* PbBr_2	White
PbF_2	White
PbI_2	Orange-yellow
PbSO_4	White
$\text{Pb}(\text{OH})_2$	White
PbC_2O_4	White
PbCO_3	White
PbCrO_4	Yellow
$\text{Pb}_3(\text{PO}_4)_2$	White
PbS	black

* Decreasing Concentration of Pb^{++} Ion

1. Chloride ion solution (use dilute HCl)

When a dilute hydrochloric acid is added to lead ions solution. a white precipitate of lead chloride, PbCl_2 is formed.




a) PbCl_2 chemical decomposition is **not affected** by ultraviolet or to a strong sun light (distinction from AgCl and Hg_2Cl_2).


b) PbCl_2 is soluble in hot water, but separates again, as white upon cooling the solution.

2. Sulfate ion solution (use dilute H_2SO_4).

Lead ion and sulfate ion combine to form a white precipitate of lead sulfate (PbSO_4)




PbSO_4 is soluble in concentrated solution of ammonium acetate.




3. Chromate ion solution (use K_2CrO_4 solution)

Lead ion and chromate ion combine to form a yellow precipitate of lead chromate, PbCrO_4 .




PbCrO_4 , which is soluble in nitric acid and in alkali hydroxides.



4. Sulfide ion solution (use H_2S which means Thioacetamide CH_3CSNH_2)

Lead ion and sulfide ion combine to form a black precipitate of lead sulfide, PbS .



PbS , which is soluble in hot dilute nitric acid.



C. Identification of the Mercurous ion Hg_2^{2+}

A list of some of the insoluble mercurous compounds is given in table 3 below:


Table 3 Insoluble mercurous compounds with their colors

Compound	color
* Hg_2Cl_2	White
$\text{Hg}_2\text{C}_2\text{O}_4$	White
Hg_2Br_2	Whitish yellow
Hg_2CO_3	Yellowish brown
Hg_2I_4	Yellow
Hg_2S	black

* Decreasing Concentration of Hg_2^{++} Ion

1. Chloride ion solution (use dilute HCl)


When a dilute hydrochloric acid is added to mercurous ions solution. a white precipitate of mercurous chloride, Hg_2Cl_2 is formed.



- When Hg_2Cl_2 exposed to ultraviolet or to a strong sun light a bright reddish orange fluorescence is noticed (distinction from PbCl_2 and AgCl).
- Hg_2Cl_2 turns to a black mixture when treated with aqueous ammonia.




This mixture of HgNH_2Cl and Hg is soluble in concentrated nitric acid.




2. Aqueous Ammonia

Mercurous ions form when react with aqueous ammonia a black precipitate of finely divided mercury



3. 8-quinolinol solution (8-hydroxy quinolone or oxime).


When this reagent is added to a solution containing mercurous ions. Gelatinous yellow precipitate of mercurous 8-quinolinolate, $\text{Hg}_2(\text{C}_9\text{H}_6\text{NO})_2$, is formed, in the present of ammonium acetate.



The precipitate changes its color from yellow to yellowish orange, to orange and then to dark orange in about five minutes.

4. Stannous chloride solution (SnCl_4 solution)

Stannous chloride reduces the mercuric ion to the mercurous ion giving a white precipitate of Hg_2Cl_2 .






Excess of the stannous chloride solution causes the reduction of mercurous chloride, first formed, to black metallic mercury).



The principle of separation

Due to the fact that only AgCl , PbCl_2 and Hg_2Cl_2 are slightly soluble in dilute acid solutions, the ions Ag^+ , Pb^{2+} and Hg_2^{2+} can be separated from the rest of the cations by the addition of chloride ions.




Lead chloride can be separated from AgCl and Hg_2Cl_2 by treating the chloride mixture with boiling water. Lead ion is then detected, in the filtrate, by the addition of the chromate ion.

Silver chloride is separated from mercurous chloride, on the basis of its solubility in aqueous ammonia which produces the complex ion, $\text{Ag}(\text{NH}_3)_2^+$. For the identification of silver ion, the solution containing this complex ion is acidified by nitric acid, thus AgCl is formed again.

Aqueous ammonia, not only dissolves AgCl but, also, converts Hg_2Cl_2 to a grayish black residue which is a mixture of HgNH_2Cl and Hg . This black residue serves as a good indication for the presence of Hg_2^{2+} .

Analysis of group I cations



The unknown solution

A 10 ml unknown solution will be provided to each student, only after the completion of the analysis of the known solution. The student should use 5 ml and save the rest to be used in case of loss of the first portion accidentally during the analysis or, when repetition of the work becomes necessary.

Procedure NO. I Separation and identification of Group I cations

<p>Precipitation of the group</p> <p>Place 10 drops of the given solution in a test tube and add 3 drops of 3M HCl (1). Heat the mixture in a water bath a few minutes, cool under the tap and allow the ppt. to settle. test for complete ppt. centrifuge and save the centrifuge , wash the ppt. with 10 drops of cold water containing a few drops of 3M HCl (2).</p>		
<p>Precipitate:</p> <p>White ppt. which may be AgCl, Hg₂Cl₂ and PbCl₂ wash twice with 7 drops portion of boiling water and combine washings (3).</p>		<p>Centrifugate Dist. water</p>
<p>Centrifugate:</p> <p>May contain Pb⁺⁺ and Cl⁻ ions. Divide, while still hot, into two parts:</p> <p>Small part:</p> <p>Add a few drops of 1M K₂CrO₄ solution, yellow ppt. of PbCrO₄ confirms presence of Pb⁺⁺ion.</p> <p>Large pare:</p> <p>Evaporate to half its volume cool to room temperature. formation of white crystalline ppt. not affected by U.V. light confirms Pb⁺⁺.</p>	<p>Residue :</p> <p>May contain AgCl and Hg₂Cl₂. Apply the U.V. light test on a small part of the ppt. (4). Blackening confirms Ag⁺ ion, but reddish orange fluorescence confirms Hg₂⁺⁺ion. If the silver ion is present the identification of Hg₂⁺⁺ion. May not be clear. In this case, treat the white ppt. with 10 drops of 3M warm aqueous ammonia, then centrifuge.</p>	
	<p>Centrifugate</p> <p>May contain Ag(NH₃)₂⁺ and Cl⁻ ions which give in nitric acid a white ppt. of AgCl indicating the presence of Ag⁺ ion.</p>	<p>Residue:</p> <p>If black, due to (Hg+HgNH₂Cl).Hg₂⁺⁺ ion is confirmed.</p>

Notes on procedure No. I

- slight excess of dilute hydrochloric acid decreases the solubility of the chloride thus, cases more complete precipitation. But large excess of a HCl or large concentration of the chloride ion should be avoided since it will partially dissolve AgCl and PbCl₂ through the formation of complex ion.

$$\begin{array}{l} \text{AgCl} + \text{Cl}^- \longrightarrow \text{AgCl}_2^- \\ \text{AgCl} + 2\text{Cl}^- \longrightarrow \text{AgCl}_2^- \\ \text{PbCl}_2 + 2\text{Cl}^- \longrightarrow \text{PbCl}_4^- \end{array}$$
- if the wash water contains no HCl, an appreciable amount of PbCl₂ will dissolving during the process of washing. The presence of hydrochloric decreases the solubility of lead chloride due to the common ion effect.
- since PbCl₂ is fairly soluble in hot water , its solution should be kept hot when the test for Pb⁺⁺ with CrO₄⁼ is performed. Otherwise lead ion may not be confirmed.
- the ultra violet test for the silver ion is very sensitive and reliable. Silver chloride is always easy to identify. But unless Hg₂Cl₂ is present in a relatively large quantity, it can not be detected since the grayish black color of the decomposed AgCl masks the reddish orange color of the fluorescent mercurous chloride.

Report sheet

Analysis of group I cations

Name:.....

Unknown No.:.....

Section:.....

Date:.....


Notes and Equations :

<i>Ions of group I</i>	<i>Ions found, by the student</i>	<i>Ions present in the unknown solution</i>	<i>grade</i>
Ag^+			
Pb^{++}			
Hg_2^{++}			


Signature of the student

Signature of the grader

Group II Cations



The cations of this group are characterized by being precipitated the sulfides from a solution which is 0.3M with respect to hydrochloric acid. The precipitating reagent for this group is either hydrogen sulfide, H_2S , or thioacetamide, CH_3CSNH_2 . The latter is a source of H_2S . It is very soluble in water and hydrolyzes rapidly in aqueous solution with the liberation of hydrogen sulfide specially at a temperature higher than room temperature.



H_2S can also be generated by gentle heating of a mixture called, Aitch-Tu-Ess which is composed of sulfur, paraffin and asbestos.

A. Identification of the Mercuric ion Hg^{++}

A list of some of the insoluble mercuric compounds is given in table 4 below:


Table 4 Insoluble mercuric compounds with their colors

Compound	color
* \downarrow $\text{Hg}(\text{IO}_3)_2$	White
HgBr_2	White
$\text{Hg}(\text{CNS})_2$	White
HgI_2	Red or yellow
\downarrow HgS	black

* Decreasing Concentration of Hg^{++} Ion

1. Sulfide ion solution (use H_2S or 1M CH_3CSNH_2 sol)

Mercuric ion and sulfide ion combine to form a black precipitate of mercuric sulfide, HgS .



HgS is soluble in sodium sulfide solution and also in aqua regia (3 volumes of concentrated hydrochloric acid and 1 volume of concentrated nitric acid).



* H_2S is written here in its molecular form instead of S^- ion, since it is slightly ionized
(i.e. weak electrolyte).


2. Stannous chloride solution

Stannous chloride reduces the mercuric ion to the mercurous ion giving a white precipitate of Hg_2Cl_2 .




(Mercuric chloride was used as the source of Hg^{++} ions, and was written in its molecular formula because it is a weak salt).

Excess of the stannous chloride solution causes the reduction of mercurous chloride, first formed, to black metallic mercury).




3. Iodide Ion solution (use KI solution).

The dropwise addition potassium iodide solution to a solution containing mercuric ions produces a red precipitate of mercuric iodide, HgI_2




HgI_2 is soluble in excess of KI solution to the formation of the mercuric-iodide complex ion, HgI_4^- .



4. Cobaltous- Thiocyanate reagent

With the mercuric ion, this reagent gives a deep blue precipitate of cobaltous mercuric thiocyanate, $\text{HgCo}(\text{CNS})_4$. The reaction is carried out by adding, to 1 ml of the solution containing Hg^{++} ions, about 1 ml of 3M sodium acetate solution and about 0.5 ml of 0.1M $\text{Co}(\text{NO}_3)_2$ solution and stirring thoroughly with dropwise addition of 0.1M ammonium thiocyanate solution.



B. Identification of the Bismuth ion Bi^{+++}

A list of some of the insoluble bismuth compounds is given in table 5 below:


Table 5 Insoluble bismuth compounds with their colors

Compound	color
* BiOCl	White
BiI_3	Reddish or brownish-gray
BiAsO_4	White
$(\text{BiO})\text{OH}$	White
BiPO_4	White
$\text{Bi}(\text{OH})_3$	White
Bi_2S_3	Brown


* Decreasing Concentration of Bi^{+++} Ion

1. Sulfide ion solution (use H₂S or CH₃CSNH₂)

Bismuth ion and sulfide ions give a brown-black precipitate of bismuth sulfide, Bi₂S₃.




Bi₂S₃ is insoluble in cold dilute acids, but soluble in boiling hydrochloric acid and in hot dilute nitric acid.



2. Aqueous ammonia




Bismuth forms a white precipitate of bismuth hydroxide, Bi(OH)₃, when reacts with aqueous ammonia.



Bi(OH)₃ is insoluble in excess of the precipitating reagent.

3. Stannite ion, SnO₂⁼, (each student prepares his solution).

Bismuth ions give a black precipitate of finely divided metallic bismuth when react with stannite ions. To prepare sodium stannite, Na₂SnO₂, solution add 3M NaOH solution drop by drop to a stannous chloride solution until the white precipitate of stannous hydroxide, Sn(OH)₂, just disappears.



C. Identification of the Cupric ion Cu⁺⁺

A list of some of the insoluble bismuth compounds is given in table 6 below:


Table 6 Insoluble bismuth compounds with their colors

Compound	color
* Cu(IO ₃) ₂	Blue
CuCrO ₄	brown
CuC ₂ O ₄	Bluish-white
CuCO ₃	Blue
Cu ₂ Fe(CN) ₆	Red-brown
Cu(OH) ₂	Blue
↓ CuS	Black


* Decreasing Concentration of Cu⁺⁺ Ion

1. Sulfide ion solution (use H₂S or CH₃CSNH₂)

Cupric ion and sulfide combine to form a black precipitate of cupric sulfide, CuS.




CuS is insoluble in sodium or potassium hydroxide, but soluble in hot dilute nitric acid.



2. Aqueous ammonia

Cupric ion forms with aqueous ammonia the deep blue color solution of tetra ammine cupric complex ion, Cu(NH₃)₄⁺⁺.




The deep blue color of the complex ion fades out upon acidifying the solution, due to the destruction of the complex ion.



3. 8-quinolinol solution (known as 8-hydroxy quinoline or oxine).


A light green precipitate of cupric 8-quinolinolate, Cu(C₉H₆NO)₂, is formed when cupric ion reacts with this reagent in a solution saturated ammonium acetate.



The precipitate has a low solubility in acetone and in ethanol.


4. Potassium ferrocyanide solution (use K₄[Fe(CN)₆] solution).

A reddish-brown precipitate is formed when cupric ions react with potassium ferrocyanide solution.



5. Iodide ion solution (use KI or NaI solution)

A yellowish- brown mixture of Cu₂I₂ (white) and I₂ (yellowish brown) precipitate is formed when cupric ions react with Iodide.



D. Identification of the Cadmium ion Cd⁺⁺

A list of some of the insoluble bismuth compounds is given in table 7 below:


Table 7 Insoluble cadmium compounds with their colors

Compound		color
* ↓	Cd(CN) ₂	White
	CdC ₂ O ₄	White
	Cd(OH) ₂	White
	Cd ₃ (AsO ₄) ₂	White
	CdCO ₃	White
	CdS	Yellow-orange

* Decreasing Concentration of Cd⁺⁺ Ion


1. Sulfide ion solution (use H₂S or CH₃CSNH₂)

Cadmium ion and sulfide combine in a neutral or a slightly acidic solution to form a yellow precipitate of cadmium sulfide, CdS. Which soluble in hot dilute nitric acid and also in hot diluted sulfuric acid but insoluble in KCN solution.



2. Aqueous ammonia

A white precipitate of cadmium hydroxide, Cd(OH)₂, is formed when cadmium ion reacts with aqueous ammonia.




Cd(OH)₂ is soluble in excess of the precipitant due to the formation of cadmium tetrammine complex ion Cd(NH₃)₄⁺⁺.



3. 8-quinolinol solution (known as 8-hydroxy quinoline or oxine).

A light greenish yellow precipitate of cadmium 8-quinolinolate, Cd(C₉H₆NO)₂.2H₂O, is formed when cadmium ion reacts with this reagent in a solution saturated ammonium acetate.



The precipitate fluoresces intensely pale green in ultra-violet light, and is insoluble in acetone or ethanol.

The principle of separation of group II cations


The copper group of the cations is composed of five cations which are precipitated as sulfides from a solution 0.2-0.3M in H_3O^+ ions using hydrochloric acid, and hydrogen sulfide or thioacetamide as the precipitating agent. Under this condition group II cations are separated from those of the subsequent groups.

Not only the sulfides of group II cations are insoluble in water but also some of the sulfides group III cations. It can be seen from the table below that all the sulfides involved can be considered as being very slightly soluble, and vary widely in their solubilities.

Insoluble sulfides of some cations from group II and group III with their solubility product constants table (8)

Sulfide from group II (insoluble in 0.3M HCl)		Sulfide from group III (soluble in 0.3M HCl)	
<u>Sulfide</u>	<u>Ks.p.</u>	<u>Sulfide</u>	<u>Ks.p.</u>
PbS	3.4×10^{-28}	MnS	1.4×10^{-15}
CdS	3.6×10^{-29}	FeS	3.7×10^{-19}
CuS	8.5×10^{-45}	CoS	7.0×10^{-23}
HgS	3.0×10^{-53}	ZnS	3.0×10^{-23}
Bi_2S_3	1.6×10^{-72}	NiS	1.4×10^{-24}
Sb_2S_3	1.6×10^{-93}		

The copper group
Analysis of group II cations



Procedure NO. 2 Separation and identification of Group II cations

<p>Precipitation of the group</p> <p>Place 10 drops of the given solution in a test tube and adjust the acidity to 0.3M HCl (1) by using 1M HCl. Add 10 drops of H₂S and then heat to boiling for a few minutes, centrifuge and test the solution for complete precipitation of the group cations as sulfides by heating and passing H₂S again stirr well and centrifuge.</p>		
<p>Precipitate:</p> <p>Add 10 drops of 3M HNO₃ upon the ppt. with stir, and heat in water bath for 3 minutes and centrifuge .</p>		<p>Centrifugate</p> <p>Groups III-V cations</p>
<p>Precipitate:</p> <p>May contain HgS (b) or Hg(NO₃)₂.2HgS (w), dissolve in 5 drops of aqua regia, add 10 drops of water, then heat for about 3 minutes (under hood). Cool then add 2 drops of SnCl₂ sol. white ppt. , Hg₂Cl₂ or gray (Hg + Hg₂Cl₂) confirms presence of Hg⁺⁺ ion.</p>	<p>Centrifugate</p> <p>Contains Pb⁺⁺, Bi⁺⁺⁺, Cu⁺⁺ and Cd⁺⁺ ions. Make the sol. basic by dropwise addition of concentrated aq. Ammonia with stirring until about 3 drops in excess of conc. NH₃ has been added (test litmus). Then centrifuge.</p>	
<p>Precipitate :</p> <p>This ppt. contains the hydroxides of Pb⁺⁺ and Bi⁺⁺⁺ as Pb(OH)₂, Bi(OH)₃, add 4 drops of NaOH and place in the boiling water for 2 min. and centrifuge.</p>		<p>Centrfugate :</p> <p>May be contains Cu(NH₃)₄⁺⁺ and Cd(NH₃)₄⁺⁺ complex ions, if solution is deep blue Cu⁺⁺ ion is confirmed (4). (If necessary concentrate the solution to 1 ml by gentle heating), divide the solution into two parts:</p>
		<p>First portion :</p> <p>Add drop by drop carefully of KCN solution until the blue color of the solution disappears (5), add 10 drops of H₂S and boil for 5 min. formation of yellow ppt. of CdS means Cd⁺⁺ ions found.</p>
<p>Centrfriugate:</p> <p>This solution contains Na₂[PbO₂] which acidify with diluted acetic acid and then add 2 drops of potassium chromate K₂CrO₄ a yellow ppt. of PbCrO₄ will appears.</p>	<p>Precipitate :</p> <p>This white ppt. contains bismuth hydroxide Bi(OH)₃, wash with 10 drops of water and prepare sodium stannite (2), (3) to sol. from which 10 drops of are poured upon the ppt. immediate blackening of the ppt. (Bi) confirms Bi⁺⁺⁺ ion.</p>	<p>Second portion:</p> <p>Acidify with diluted acetic acid and add 1 drop of potassium ferrocyanide solution, K₄[Fe(CN)₆]. The formation of reddish-brown precipitate means the Cu⁺⁺ ions is found.</p>

Notes on procedure No. II

- 1- To adjust the acidity to about 0.3M HCl, add 1 drop of Methyl Violate dye to the mixture solution, the Methyl Violate shows different colors at different hydronium ion concentration (different acidity) as shown below:

Color	Solution
Violet	neutral and basic
Blue	0.1 M H_3O^+
Blue-green	0.3 M H_3O^+
Yellow-green	0.5 M H_3O^+
Yellow	1.0 M H_3O^+ and higher

Blue $\xrightarrow[1 \text{ drop HCl}]{\text{Green (0.3M)}}$ $\xleftarrow[1 \text{ drop NH}_4\text{OH}]{\text{Yellow}}$


If the color of solution is blue,
The solution will be alkaline

If the color of solution is yellow,
The solution will be acidic

- 2- to prepare sodium stannite solution (see the preparation on the Identification of the Bismuth ion Bi^{+++} no. 3) this reagent must be fresh, when used, and made up just before using it, since it undergoes an atmospheric oxidation to sodium stannate, Na_2SnO_3 , which is not oxidized by $\text{Bi}(\text{OH})_3$.
- 3- the stannite ion, SnO_2^- , reduces bismuth hydroxide to metallic bismuth, Bi , which is black. Thus bismuth ion is not confirmed unless the precipitate is black and not brown.
- 4- If the blue color is not quite satisfactory, place 1 ml of the solution in a test tube then acidify with 3M HCl and add a solution of $\text{K}_4[\text{Fe}(\text{CN})_6]$. the formation of reddish-brown precipitate of $\text{Cu}_2[\text{Fe}(\text{CN})_6]$ means the Cu^{++} ion is found.
- 5- potassium cyanide, KCN , solution converts the cupric-ammonia complex ion, $\text{Cu}(\text{NH}_3)_4^{++}$, to the stable and colorless complex ion, $\text{Cu}(\text{CN})_3^-$, the later dose not react with hydrogen sulfide, thus no black precipitate of CuS is formed.

Systematic group separation of the cations

The overall scheme for the systematic analysis of a mixture of cations, based on the five major groups discussed above, is shown as a flow chart below.



Report sheet

Analysis of group II cations

Name:.....

Unknown No.:.....

Section:.....

Date:.....

Notes and Equations :

<i>Ions of group II</i>	<i>Ions found, by the student</i>	<i>Ions present in the unknown solution</i>	<i>grade</i>
Pb ⁺⁺			
Hg ⁺⁺			
Bi ⁺⁺⁺			
Cu ⁺⁺			
Cd ⁺⁺			

Signature of the student

Signature of the grader

Group III
The Aluminum-Zinc Group
Al⁺⁺⁺, Cr⁺⁺⁺, Fe⁺⁺⁺
Fe⁺⁺, Co⁺⁺, Ni⁺⁺, Mn⁺⁺ and Zn⁺⁺

This group is composed of eight cations whose chlorides and sulfides are soluble in acidic solution. These cations, above are separated from group IV and V as a mixture of insoluble sulfides and hydroxides in a solution containing ammonia and ammonium chloride, using NH₃ and H₂S or (NH₄)₂S as the precipitating agents.

The insoluble sulfides of group III cations have relatively high solubility product constants, compared with those of group II, (see the table 8) that prevent their precipitation with group II, and the hydroxides have smaller solubility product constants than their corresponding sulfides. Therefore, two precipitating reagents are used for the precipitation and separation of cations of this group from those groups IV and V. ammonium sulfide or Fe⁺⁺, Co⁺⁺, Ni⁺⁺, Mn⁺⁺ and Zn⁺⁺ ions, while aqueous ammonia precipitates, as hydroxides the cations Al⁺⁺⁺, Cr⁺⁺⁺, Fe⁺⁺⁺

A . Identification of the aluminum ion, Al⁺⁺⁺

A list of some of the insoluble aluminum compounds is given in table 9

Table 9 Insoluble aluminum compounds with their colors



Compound	color
* ↓ Al ₂ S ₃	White
Al(OH) ₃	White
AlAsO ₄	White
↓ AlPO ₄	White

* Decreasing Concentration of Al⁺⁺⁺ Ion

1. Hydroxide ion solution

a. Aqueous ammonia


A white precipitate of aluminum hydroxide, Al(OH)₃, is formed when aluminum ions react with aqueous ammonia. The precipitate is rather Al₂O₃.xH₂O



The precipitate is insoluble in excess of the precipitating reagent, but soluble in acidic media.


b. Alkali hydroxide solutions.

Aluminum hydroxide is also formed as above but excess of the alkali hydroxide solutions dissolves the precipitate.



2. 8-quinolinol solution (known as 8-hydroxy quinoline or oxine).

a yellowish green precipitate of 8-quinolinolate, $\text{Al}(\text{C}_9\text{H}_6\text{NO})_3$, is formed in a solution saturated ammonium acetate .



3. Aluminon solution

The addition of aluminon solution to a basic solution with ammonia containing aluminum ions and ammonium acetate, causes the adsorption of this red organic dye (aluminon) by aluminum hydroxide. No precipitate is actually formed, but a jelly type red precipitate (lake) is produced.

B. Identification of the chromium ion, Cr^{+++}

A list of some of the insoluble chromium compounds is given in table 10

Table 10 Insoluble chromium compounds with their colors


Compound	color
* $\text{Cr}(\text{OH})_2$	Yellowish-brown
$\text{Cr}(\text{OH})_3$	Yellowish-green
CrAsO_4	Light green
CrPO_4	Green With (4H ₂ O) or violate (3 OR 6H ₂ O)

* Decreasing Concentration of Cr^{+++} Ion

1. Hydroxide ion solution

a. Aqueous ammonia

A gelatinous grayish-green precipitate of chromic hydroxide, $\text{Cr}(\text{OH})_3$, is formed when chromic ion react with aqueous ammonia. The precipitate is rather $\text{Cr}_2\text{O}_3 \cdot x\text{H}_2\text{O}$




$\text{Cr}(\text{OH})_3$ is readily dissolves in acids, but slightly soluble in excess aqueous NH_3 in the cold giving a violate or pink solution of the chromic complex ion, $\text{Cr}(\text{NH}_3)_6^{+++}$

$$\text{Cr}(\text{OH})_3 + 6\text{NH}_3 \text{ aq.} \rightleftharpoons \text{Cr}(\text{NH}_3)_6^{3+} + 3\text{OH}^-$$


b. Alkali hydroxide solutions.

A precipitate similar to the above is formed which is readily soluble in acids and in excess of the alkali hydroxide solutions.




2. Sodium peroxide Na_2O_2 , (or using H_2O_2 and NaOH sol.)


Chromium ion, in the presence of sodium peroxide, is oxidized to the yellow chromate ion $\text{CrO}_4^{=}$.




The chromate ion is precipitated as a yellow PbCrO_4 or BiCrO_4



The addition of H_2O_2 to an acid solution of the chromate ion results in the productions of a blue coloration owing to the formation of CrO_5 .



CrO_5 is unstable in aqueous solution, and readily with the liberation of oxygen.



The test for chromium ion may be done as follows:

Precipitate the $\text{CrO}_4^{=}$ ion as PbCrO_4 . dissolve the precipitate in hot 3M HNO_3 , cool to room temperature and add equal volume of amyl alcohol followed by about 1 ml 3% H_2O_2 . Shake twice and notice the blue color in the alcohol layer which will fade gradually.

C. Identification of the Ferric ion , Fe^{+++}


A list of some of the insoluble ferric compounds is given in table 11 below:

Table 11 Insoluble ferric compounds with their colors


Compound		color
*	$\text{Fe}_2(\text{SeO}_3)_3$	Deep blue Reddish brown Yellowish white black
	$\text{Fe}_4[\text{Fe}(\text{CN})_6]_3$	
	Fe(OH)_3	
	FePO_4	
	↓ Fe_2S_3	

1. Hydroxide ion solution (NH₃ aq. or alkali hydroxide sol.)

Ferric ion combines with the hydroxide ions from a reddish-brown precipitate of ferric hydroxide, Fe(OH)₃, which is insoluble in excess of the precipitant, but readily dissolves in acids.




The reddish-brown precipitate does not, actually, correspond to the formula, Fe(OH)₃ but rather to the formula Fe₂O₃.xH₂O as it is the case with Al₂O₃.xH₂O and Cr₂O₃. xH₂O.




2. Sulfide ion solution (H₂S or (NH₄)₂S)

In the presence of hydrochloric acid, hydrogen sulfide reduces ferric ions to ferrous ions with the formation of free sulfur.



In a basic solution, hydrogen sulfide or ammonium sulfide solution gives a black precipitate of ferric sulfide, Fe₂S₃




The precipitate is soluble in dilute hydrochloric acid, which results in the reduction of iron from the ferric state to the ferrous state.




3. Thiocyanate ion ,CNS⁻, solution (use NH₄CNS or KCNS sol.)

Ferric ion reacts with excess thiocyanate ions to give a deep reddish brown complex ion, Fe(CNS)₆³⁻ or Fe CNS⁺⁺.



4. 8-quinolinol solution.

Fine dark green precipitate of ferric 8-quinolinolate, Fe(C₉H₆NO)₃, is formed upon adding this reagent to a solution containing ferric ions and saturated with ammonium acetate.



Ferric ion the only cation in the precipitate entire five groups of the cations that gives such color or the 8-quinolinolate.

D. Identification of the Ferrous ion, Fe⁺⁺

A list of some of the insoluble ferrous compounds is given in table 12 below:

Table 12 Insoluble ferrous compounds with their colors


Compound	color
* FeC ₂ O ₄	Yellow
Fe(OH) ₂	Pale green or white
FeCO ₃	Gray
↓ FeS	Brown-black

* Decreasing Concentration of Fe⁺⁺ Ion


1. Hydroxide ion solution

a. Aqueous ammonia

In the complete absence of air. Ferrous ion and hydroxide ions combine to form a white precipitate of ferrous hydroxide, Fe(OH)₂




However, the precipitate may look reddish-brown or muddy green, instead of white, due to the quick oxidation of Fe(OH)₂ to Fe(OH)₃ by the air oxygen.



The precipitate of Fe(OH)₂, is insoluble in excess of the precipitate, but soluble in acids.


b. Alkali hydroxide solutions.

Similar results are obtained as in above.



2. Sulfide ion solution (H₂S or (NH₄)₂S)

No precipitate is formed in acidic solutions when ferrous and sulfide ions are brought together. But a black precipitate of ferrous sulfide, FeS, is obtained from basic solutions. The precipitate turns brown upon exposure to air.




The precipitate is soluble in acids giving off hydrogen sulfide.

3. Ferrocyanide ion, Fe(CN)₆³⁻ solution (use K₄Fe(CN)₆ sol.)


If the reaction between this reagent and the ferrous ion is carried out in complete absence of air, a white precipitate of potassium ferrous ferrocyanide, K₂[Fe(Fe(CN)₆], is formed. However, partial oxidation by atmospheric oxygen may take place, which in turn produces a light blue mixture of potassium ferric ferrocyanide KFe(Fe(CN)₆].



If potassium ferricyanide, $\text{K}_3[\text{Fe}(\text{CN})_6]$, solution is used a dark blue precipitate is obtained.



In a basic solution, hydrogen sulfide or ammonium sulfide solution gives a black precipitate of ferric sulfide, Fe_2S_3



The blue precipitate turns reddish brown if sodium or potassium hydroxide solution is added to it, causing the decomposition of the precipitate and the formation of $\text{Fe}(\text{OH})_3$.

E. Identification of the Cobaltous ion, Co^{++}

A list of some of the insoluble cobaltous compounds is given in table 13 below:


Table 13 Insoluble cobaltous compounds with their colors

Compound		color
* ↓	CoC_2O_4	Raddish white
	$\text{Co}(\text{OH})_2$	Rose-red
	CoCO_3	Red
	CoS	Black

* Decreasing Concentration of Co^{++} Ion

1. Aqueous ammonia

A blue basic precipitate, $\text{Co}(\text{OH})\text{NO}_3$ is formed when aqueous ammonia is added to a solution of cobaltous nitrate.





$\text{Co}(\text{OH})\text{NO}_3$ is readily dissolves in excess of aqueous ammonia, due to the formation of the cobaltose ammonia complex ion, $\text{Co}(\text{NH}_3)_6^{++}$.




2. Alkali hydroxide solutions.

A pink precipitate of cobaltous hydroxide, $\text{Co}(\text{OH})_2$ is formed in the presence of slight excess NaOH solution. $\text{Co}(\text{OH})_2$ is soluble in excess of the alkali hydroxide, if the concentration of the hydroxide ions is higher than 8 molar.




Exposure of Co(OH)_2 to atmospheric oxygen leads to the formation of brown cabaltic hydroxide, Co(OH)_3 .

Addition of H_2O_2 or Na_2O_2 in a high concentration of the hydroxide ions converts Co(OH)_2 to Co(OH)_3 very rapidly.




3. Sulfide ion solution (H_2S or $(\text{NH}_4)_2\text{S}$)

In a neutral or basic solution with ammonia, the sulfide ion combines with the cobaltuse ion to form a black precipitate of cobaltose sulfide, CoS .




Cobaltose sulfide in insoluble in a very dilute acids and in acetic acid, but dissolved readily by concentration HCl , HNO_3 and by aqua regia. In the case of HCl , CoS dissolves quite slowly.



4. Thiocyanate ion

when an alcoholic solution of NH_4CNS is added slowly to the surface of a slightly acidic solution containing Co^{++} ions without mixing, an intense blue color is formed, because of the formation of the cobalti thiocyanate complex ion, Co(CNS)_4^- the color appears at the junction of the two liquids.



Mixing or diluting the solution, containing Co(CNS)_4^- ions, with water destroys this complex ion.

The same test may be performed in a different way, by adding concentrated NH_4CNS solution then extracting the blue complex ions of Co(CNS)_4^- by shaking with amyl alcohol.

F. Identification of the Nickelous ion, Ni^{++}

A list of some of the insoluble nickelous compounds is given in table 14 below:


Table 14 Insoluble nickelous compounds with their colors

Compound		color
* ↓	$\text{Ni(IO}_3)_2$	Yellow
	NiCO_3	Light green
	Ni(OH)_2	Green
	Ni(CN)_2	Green
	NiS	Black


* Decreasing Concentration of Ni^{++} Ion

1. Alkali hydroxide solutions.

A solution of NaOH . For example, precipitates the light green gelatinous nickelous hydroxide, Ni(OH)_2 from solutions of Ni^{++} ion.




Ni(OH)₂ is soluble in solutions containing NH₃ or NH₄⁺ ions producing a greenish blue solutions of the complex ions, Ni(NH₃)₄⁺⁺ and Ni(NH₃)₆⁺⁺




Ni(NH₃)₆⁺⁺ is more stable (K_{inst}= 1.8 x 10⁻⁹) than Ni(NH₃)₄⁺⁺ (K_{inst}= 1.0 x 10⁻⁸).

2. Aqueous ammonia

A green basic precipitate is formed upon the addition of aqueous ammonia to a solution of nickelous chloride.




NiOHCl is soluble in excess of aqueous ammonia, due to the formation of the complex ion, Ni(NH₃)₆⁺⁺.




3. Sulfide ion solution (H₂S or (NH₄)₂S)

In basic solutions, the sulfide ion combines with the nickelous ion to form a black precipitate of nickelous sulfide, NiS.




The hydronium ions formed are removed from the solution through the reaction with base NH₃ in the solution, giving NH₄⁺ ions. In other words the reaction may be written as follows:-



4. Dimethylglyoxime, C₄H₈O₂N₂ solution

When an alcoholic solution of dimethylglyoxime is added to a slightly ammoniacal solution (basic with ammonia) containing Ni⁺⁺ ions, a red precipitate of nickel dimethylglyoxime, Ni(C₄H₇O₂N₂)₂ is formed.



G. Identification of the Manganous ion, Mn⁺⁺

A list of some of the insoluble Manganous compounds is given in table 15 below:

Table 15 Insoluble manganous compounds with their colors


Compound	color
Mn(OH) ₂	White pink
MnCO ₃	Rose-pink
* MnS	Green or pink
↓ MnSO ₄	white

* Decreasing Concentration of Mn⁺⁺ Ion

1. Hydroxide ion solution

a. Aqueous ammonia


Manganous ions are partially precipitated by aqueous ammonia, as white manganous hydroxide, $\text{Mn}(\text{OH})_2$.




The precipitate is soluble in solutions containing as appreciable amount of NH_4^+ ions, since the latter upon hydrolysis furnishes H_3O^+ ions which will neutralize the OH^- ions, thereby reducing the OH^- ions in solution which are in equilibrium with $\text{Mn}(\text{OH})_2$.

b. Alkali hydroxide solutions .




Addition of NaOH solution, for example, to a solution containing Mn^{++} ions results in the complete precipitating of white manganous hydroxide, $\text{Mn}(\text{OH})_2$



The precipitate is insoluble in excess of the precipitant. Exposure of $\text{Mn}(\text{OH})_2$ to the air results in the conversion of the precipitate to the brown manganic hydroxide, $\text{Mn}(\text{OH})_3$ or the hydrated dioxide, $\text{MnO}_2 \cdot x\text{H}_2\text{O}$.




Manganous ion, in the presence of Na_2O_2 and excess NaOH produces a brown precipitate of the hydrated manganese dioxide.



3. Sulfide ion solution (H_2S or $(\text{NH}_4)_2\text{S}$ sol.)


In a slightly basic solutions, manganous ion gives a light pink precipitate of, manganous sulfide when reacts with the sulfide ion.



MnS is soluble in dilute acids and in acetic acid this is because MnS is not very insoluble, thus it furnishes a relatively high concentration of the sulfide ions that are removed by the formation of HS^- ion and H_2S .

4. Bismuthate ion solution (use NaBiO_3).

Manganous ions, in a cold solution containing dilute nitric acid are oxidized by the bismuthate ions (usually added as solid NaBiO_3) to the reddish purple or violet-red permanganate ions MnO_4^- .



High concentration of the Mn^{++} ion may produce dark solution, but the color can be clearly noticed after filtering off the excess of solid bismuthate and then diluting the filtrate with water.

Several other oxidizing agents may be used for the oxidation of Mn^{++} ion to MnO_4^- ion, such as PbO_2 , $(\text{NH}_4)_2\text{S}_2\text{O}_8$ and KIO_4 but with different procedures.

H. Identification of the Zinc ion, Zn^{++}

A list of some of the insoluble zinc compounds is given in table 16 below:

Table 16 Insoluble zinc compounds with their colors


Compound		color
* ↓	$\text{Zn}(\text{CN})_2$	White
	ZnC_2O_4	White
	ZnSO_3	White
	$\text{Zn}_2\text{Fe}(\text{CN})_6$	White
	$\text{Zn}(\text{OH})_2$	White
	ZnS	White

* Decreasing Concentration of Zn^{++} Ion

1. Hydroxide ion solution


a. Aqueous ammonia

A white precipitate of zinc hydroxide, $\text{Zn}(\text{OH})_2$ is formed upon the addition of aqueous ammonia, to a solution of zinc ion.




b. Alkali hydroxide solutions .


Zinc ion reacts with the alkali solution giving a white gelatinous precipitate of zinc hydroxide, $\text{Zn}(\text{OH})_2$



$\text{Zn}(\text{OH})_2$ is soluble in excess of the precipitant,




And also soluble in dilute acids,




3. Sulfide ion solution (H_2S or $(\text{NH}_4)_2\text{S}$)


In a basic solutions, zinc ion gives a white precipitate of, zinc sulfide (ZnS) when reacts with the sulfide ion in aqueous ammonia solution.



Partial precipitation, as sulfide, take place in neutral or very slightly acidic solutions. If the ammoniacal solution of Zn^{++} ions is acidified with acetic acid can then saturated with H_2S , the precipitate of ZnS will be formed.



ZnS is insoluble in excess of the precipitant, in acetic acid and in alkali solutions, but readily dissolves in mineral acids.



4. 8-quinolinol solution.

Zinc ion and 8-quinolinolate react with each other, in a solution saturated with ammonium acetate to give a light yellow precipitate of zinc 8-quinolinolate $\text{Zn}(\text{C}_9\text{H}_6\text{NO})_2$.



The principle of the separation of group III

It was mentioned that the chlorides and sulfides of group III cations are soluble in acid media, otherwise they would precipitate with either the members of group I or II or may be with both.


The separation of group III cations from groups IV and V is based upon the fact that the sulfides and hydroxides of this group are insoluble in moderately ammoniacal solution, while those of group IV and V are soluble.

Fairly complete precipitation of the members of this group is usually made in an ammoniacal solution buffered with NH_4Cl using both the sulfide and hydroxide ions as the precipitants. The presence of NH_3 in the solution has two functions. First, it will precipitate the hydroxides (or rather the hydrous oxides) of Al, Cr and Fe, and second it will increase the concentration of the sulfide ions to a value that permits the solubility product constants of ZnS, FeS, NiS, CoS and MnS to be exceeded.


The hydroxide ion concentration must be sufficiently high to precipitate the above hydroxides, and to increase the sulfide ion concentration so that the above cations are precipitated. But it should not be so high as to cause the precipitation of magnesium ion (a member of group IV cations) as magnesium hydroxide, $\text{Mg}(\text{OH})_2$.

The insoluble hydroxides and sulfides of group III cations are soluble in concentrated hydrochloric acid, with exception of CoS and NiS which are soluble in concentrated HNO_3 . When dissolution is complete, the group can then be divided into two subgroups, by the action of excess KOH or NaOH solution and Na_2O_2 . Subgroup IIIA will include aluminum, chromium and zinc ions. At the step of separation into the two subgroups, these three cations will be present as the aluminate, AlO_2^- , zincate, ZnO_2^- , and chromate, CrO_4^- ions whereas the ferric, manganous, cobaltic and nickelous ions are precipitated as their corresponding hydroxides to form subgroup IIIB.

Group III-IV and V cations




Analysis of subgroup III A cations



Analysis of subgroup III B cations

Precipitate of subgroup III B cations

$\text{MnO}_2 \cdot x\text{H}_2\text{O}$, $\text{Fe}_2\text{O}_3 \cdot x\text{H}_2\text{O}$, $\text{Co}(\text{OH})_3$, $\text{Ni}(\text{OH})_2$
or $\text{Mn}(\text{OH})_2$, $\text{Fe}(\text{OH})_3$, $\text{Co}(\text{OH})_3$ and $\text{Ni}(\text{OH})_3$



* DMG = dimethylglyoxime = $\text{C}_4\text{H}_8\text{O}_2\text{N}_2$

Procedure NO. 3 Analysis of Group III cations

Precipitation of the group:

The centrifugate (from group II, procedure No.2), may contain any or all the cations of groups III. IV and V. adjust the volume of the centrifugate 1 ml. Add 4 drops of saturated NH_4Cl solution then concentrated NH_3 (1) solution drop by drop until the solution is just basic. Add 1 ml of the bases and pass H_2S until precipitation is complete centrifuge.

Precipitate:

Contains FeS (b), CoS (b), NiS (b), MnS (p), ZnS (w), $\text{Al}(\text{OH})_3$ (w), $\text{Fe}(\text{OH})_3$ (r-br), $\text{Cr}(\text{OH})_3$ (gr-grn). Wash twice with water and discard washing. Treat the ppt. with 5 drop of concentrated HCl .

Stir thoroughly for about a minute. If the ppt. dissolves completely, Ni and Co must be absent (3). Add 3 drops of concentrated HNO_3 and heat until the precipitate dissolves. Make the solution basic with 3M KOH sol. then , slowly and carefully, add 0.2 gram of Na_2O_2 (4), (5). Heat for a few minutes and centrifuge. Wash the residue with 10 drops of water and combine washing with centrifuge.

Centrifugate

Contains groups IV and V Acidify immediately with hydrochloric acid and boil (2) to 1 ml , discarding any ppt. which may form reserve the colorless sol. for the analysis of groups IV and V cations.

Precipitate of subgroup IIIB cations



Contains $\text{MnO}_2 \cdot x\text{H}_2\text{O}$, $\text{Co}(\text{OH})_3$, $\text{Ni}(\text{OH})_2$ and $\text{Fe}(\text{OH})_3$. Wash with 10 drops of water treat according to procedure No. 2

Centrifugate (subgroup IIIA cations)

Contains AlO_2^- , ZnO_2^- , and CrO_4^- ions. Treat according to procedure No.1

Notes on procedure :

1. Aqueous ammonia precipitates some of the cations as hydroxides or oxides, and also neutralizes the hydronium ions which are used by the reaction of the cations with H_2S . for example:



In the solution is not kept basic, complete precipitation may be prevented.

Ammonia and ammonium chloride from a buffer solution which represses the OH^- ion concentration to certain values at which $\text{Mg}(\text{OH})_2$ cannot be precipitated.

2. Unless the sulfide ion is removed from the solution by acidifying and boiling ,



It will be oxidized to the sulfate ion by the atmospheric oxygen.



The formed $\text{SO}_4^{=}$ ion may cause a loss of some cations of group IV such as Sr^{++} and Ba^{++} because of the formation of SrSO_4 and BaSO_4 which are both insoluble.

3. Concentrated HCl , readily dissolves $\text{Al}(\text{OH})_3$, $\text{Cr}(\text{OH})_3$, $\text{Fe}(\text{OH})_3$, FeS , MnS and ZnS but very slowly the precipitate of CoS and NiS . A black ppt. after the addition of concentrated hydrochloric acid which dissolves in concentrated HNO_3 may indicate any or both, CoS and NiS .

4. A large excess of Na_2O_2 in a strongly alkaline solution may partially disperse $\text{Fe}(\text{OH})_3$, $\text{Co}(\text{OH})_3$ and $\text{Ni}(\text{OH})_2$ which in turn affect the test for Al^{+++} and Zn^{+++} ions.

5. sodium peroxide, Na_2O_2 , oxidizes the chromic ion, Cr^{++} , in the presence of KOH or NaOH to the chromate ion, $\text{CrO}_4^{=}$, $\text{Fe}(\text{OH})_2$ to $\text{Fe}(\text{OH})_3$ and $\text{Co}(\text{OH})_2$ to $\text{Co}(\text{OH})_3$.

Procedure NO. 1 Analysis of subgroup III A cations

Centrifugate, (from procedure of subgroup IIIA cations)

The centrifugate contain AlO_2^- , $\text{ZnO}_2^{=}$, and $\text{CrO}_4^{=}$ ions. Add 10 drops of saturated ammonium acetate solu. Then add drop by drop of acetic acid solution until the test solution becomes acidity and then add aqueous ammonia until the solution becomes ammoniacal. Then centrifugate.

Precipitate:

May contain gelatinous white ppt. $\text{Al}(\text{OH})_3$, wash with 10 drops water, then add 3M HCl drop by drop until the ppt. disappeared, add 3 drops of saturated

Centrifugate

May contain $\text{Zn}(\text{NH}_3)_4^{++}$ and $\text{CrO}_4^{=}$ add dropwise 1M of barium acetate $\text{Ba}(\text{C}_2\text{H}_3\text{O}_2)_2$ solution until precipitation is complete. Then centrifuge.

<p>ammonium acetate solution and divided into two unequal parts.</p> <p>Large porting: Add 2-3 drops of aluminon solution then render basic with 3M of aq. ammonia. Formation of red ppt. confirms presence of Al^{+++} ion.</p> <p>Smaller part: Make the solution almost neutral with 3M aq. ammonia, then add 5 drops of 8-quinolinol solution. Formation of a greenish yellow ppt. which is little soluble in acetone or ethanol and highly fluorescent, confirms presence of Al^{+++} ion.</p>	<p>Precipitate, $BaSO_4$ (1) $BaCrO_4$ (y), add 10 drops of HCl and warm for 1 min. in a Water bath. Centrifuge and discard the white ppt. ($BaSO_4$). Add 10 drops of amyl alcohol and 5 drops of H_2O_2 and shake well. A blue upper layer (2) containing CrO_5 confirms presence of Cr^{+++} ion.</p>	<p>Centrifugate, Contains $Zn(NH_3)_4^{++}$ acidify with 3M CH_3COOH and divided into two unequal parts: Large porting: Saturate with H_2S and warm in a water bath for few minutes, a white ppt. of ZnS confirms presence of Zn^{++} ion. If the ppt. is not white (3) or fails to form, use then the smaller part. Smaller part: Add 3 drops of saturated ammonium acetate solution then add 4 drops of 8-quinolinol solution. Yellow ppt. of $Zn(C_9H_6NO)_2$ insoluble in acetone or ethanol it confirms Zn^{++} ion.</p>
---	--	--

Notes on procedure :

1. Na_2O_2 oxidizes any $S^{=}$ ions in the solutions to $SO_4^{=}$ ions. The sulfate ions will precipitate as $BaSO_4$ when barium acetate solution is added.
2. the blue coloration, due to CrO_5 , may not last long since CrO_5 is unstable , decomposing with the evolution of oxygen gas.
3. the precipitate of ZnS may not be white but brown or black due to contamination with CoS or PbS . Lead sulfide, PbS , is formed here after some of the lead ions have escaped from group II cations.

Procedure NO. 2 Analysis of subgroup IIIB cations

<p>Precipitate:, (from procedure of subgroup IIIB cations) May contain $MnO_2 \cdot xH_2O$, $Ni(OH)_2$, $Co(OH)_3$ and $Fe(OH)_3$. wash with 10 drops of water and discard washings. Dissolve in concentrated hydrochloric acid then dilute with water to about 2ml divided equally into four portions.</p>			
1 st. portion Test for Fe^{+++}	2 nd. Portion Test for Mn^{++}	3ed. Portion Test for Co^{++}	4 th. Portion Test for Ni^{++}

<p>Divided into two parts a and b</p> <p>Part a.</p> <p>Add 2 drops of 3M NH_4CNS sol. deep red coloration proves presence of Fe^{+++}.</p> <p>Part b.</p> <p>Make solution almost neutral with 3M aqueous ammonia, then add 3 drops of 8-quinolinol solution. Deep greenish-black coloration, or fine precipitate (1) proves presence of Fe^{+++} ion</p>	<p>Add 0.2 gram of NaBiO_3 and shake well for a minute, then centrifuge, reddish-purple centrifugate, MnO_4^- ion proves presence of Mn^{++} ion.</p>	<p>Add with stirring enough solid NaF (2) that undissolved amount remains. Carefully let 10 drops alcoholic NH_4CNS sol(3), rundown side of the test tube containing the test solution. Formation of a blue-green ring proves presence of Co^{++} ion.</p>	<p>Add with stirring enough solid NaF (2) that undissolved amount remains. Add 3M aq. ammonia drop by drop until the solution is slightly basic. Add 4 drops of dimethylglyoxime sol. formation of a red precipitate proves presence of Ni^{++} ion (4).</p>
---	---	---	--

Notes on procedure :

1. if only a greenish-black coloration appears, one can get the precipitate by saturating the solution with ammonium acetate.
2. Any amount of ferric ions must be removed, since Fe^{+++} ions interfere with the test for cobalt and nickel ions. The addition of NaF results in the formation of the stable and colorless ion, FeF_6^{---} which does not react with either CNS^- ions or dimethylglyoxime.
3. The blue-cobaltothiocyanate complex ion, $\text{Co}(\text{CNS})_4^-$, is stable in alcoholic solution but not in aqueous solution. It may be destroyed during the addition of alcoholic thiocyanate solution, if the test solution is shaken or stirred.
4. in the presence of cobalt ions, which also react with dimethylglyoxime to form a dark colored complex ions, one should add an excess of the precipitant to ensure the formation of nickel dimethylglyoximate, $\text{Ni}(\text{C}_4\text{H}_7\text{N}_2\text{O}_2)_2$.

Report sheet

Analysis of group III cations

Name:.....

Unknown No.:.....

Section:.....

Date:.....

Notes and Equations :

<i>Ions of group III</i>	<i>Ions found, by the student</i>	<i>Ions present in the unknown solution</i>	<i>grade</i>
Al⁺⁺⁺			
Cr⁺⁺⁺			
Fe⁺⁺⁺			
Co⁺⁺			
Ni⁺⁺			
Mn⁺⁺			
Zn⁺⁺			

Signature of the student

Signature of the grader

Group IV
The Alkaline Earth Group
Ba⁺⁺, Sr⁺⁺, Ca⁺⁺ and Mg⁺⁺ ions

This group is composed of the ions of the alkaline earth metals and magnesium. The chlorides, sulfides and hydroxides of those cations are soluble under the conditions of the precipitations of the first three groups of the cations. These ions are characterized by the formation of a relatively large number of insoluble salts, with various anions, but the only anion which will precipitate these four cations satisfactorily is the phosphate anion.

Due to the fact that the metals of these cations are included in the same group of the periodic table, they have the same valence and resemble each other closely in their chemical reactions. These cations are also known to be excellent reducing agents, do not hydrolyze and show little tendency toward complex formation.

A . Identification of the Barium ion, Ba⁺⁺

A list of some of the insoluble barium compounds is given in table 17


Table 17 Insoluble barium compounds with their colors

Compound	color
BaF ₂	White
* BaSO ₃	White
BaC ₂ O ₄	White
↓ Ba(IO ₃) ₂	White
BaCO ₃	White
BaCrO ₄	Yellow
BaSO ₄	White
Ba ₃ (PO ₄) ₂	White


* Decreasing Concentration of Ba⁺⁺ Ion

1. Phosphate ion solution (use (NH₄)₂HPO₄).

Diammonium acid phosphate solution precipitates barium ions, from ammoniacal solution as a white gelatinous precipitate of barium phosphate, Ba₃(PO₄)₂.




Ba₃(PO₄)₂ is soluble in mineral acids, and in the absence of large excess of the acetate ions, barium phosphate is also soluble on acetic acid.



2. Sulfate ion solution (use H₂SO₄)

The sulfate ion precipitates the barium ion as a white finely divided precipitate of barium sulfate, BaSO₄.




BaSO₄ is insoluble in dilute acids, but can be dissolved by boiling with concentrated sulfuric acid where the bisulfate ion, HSO₄⁻ is formed.




3. Chromate ion solution (use K₂CrO₄ or Na₂CrO₄ sol.)

The chromate ion precipitates the barium ion as a yellow precipitate of barium chromate, BaCrO₄.



BaCrO₄ is soluble in mineral acids, but insoluble in acetic acid (distinction from SrCrO₄ and CaCrO₄).



4. Flame test:

Barium salts, especially the chloride, impart a yellow green color to the flame when heated in Bunsen flame. Clean the platinum wire by dipping it into concentrated hydrochloric acid and then heating it in the hottest oxidizing of the flame. Repeat this process until any color disappears the loop of the Pt wire in the solution being tested and then heat in the oxidizing zone of the flame. A yellowish-green coloration of the flame confirms presence of barium ion.

B . Identification of the Strontium ion, Sr⁺⁺

A list of some of the insoluble strontium compounds is given in table 18


Table 18 Insoluble strontium compounds with their colors

Compound	color
* Sr(OH) ₂	White
SrCrO ₄	Yellow
SrF ₂	White
SrSO ₄	White
SrC ₂ O ₄	White
SrCO ₃	White
↓ Sr ₃ (PO ₄) ₂	White


* Decreasing Concentration of Sr⁺⁺ Ion

1. Phosphate ion solution (use $(\text{NH}_4)_2\text{HPO}_4$).

Diammonium acid phosphate solution precipitates strontium ions, from ammoniacal solution as a white gelatinous precipitate of strontium phosphate, $\text{Sr}_3(\text{PO}_4)_2$.




$\text{Sr}_3(\text{PO}_4)_2$ is soluble in mineral acids, and in the absence of large excess of the acetate ions are absent.



2. Sulfate ion solution (use H_2SO_4)


The sulfate ion reacts with the strontium ion to product a white precipitate of strontium sulfate, SrSO_4 .




SrSO_4 slightly soluble in boiling hydrochloric acids, but sufficiently soluble to give a flame test for the strontium ion.

3. Carbonate ion solution (use $(\text{NH}_4)_2\text{CO}_3$ or Na_2CO_3 sol.)

The carbonate ion precipitates from slightly basic solutions the strontium ion as a white precipitate of strontium carbonate, SrCO_3 .



SrCO_3 is soluble in mineral acids, acetic acid and in water saturated with carbene dioxide.



4. Flame test:

Strontium salts, especially the chloride, when heated in Bunsen flame. Import a characteristic crimson (brilliant red) color to the flame.

C . Identification of the Calcium ion, Ca^{++}

A list of some of the insoluble Calcium compounds is given in table 19


Table 19 Insoluble calcium compounds with their colors

Compound	color
CaCrO_4	Yellow
$\text{Ca}(\text{OH})_2$	White
* CaSO_4	White
CaF_2	White
$\text{CaSO}_3 \cdot 2\text{H}_2\text{O}$	White
CaCO_3	White
CaC_2O_4	White
$\text{Ca}_3(\text{PO}_4)_2$	White


* Decreasing Concentration of Ca^{++} Ion

1. Phosphate ion solution (use $(\text{NH}_4)_2\text{HPO}_4$).

Diammonium acid phosphate solution precipitates calcium ions, from ammoniacal solution as a white gelatinous precipitate of calcium phosphate, $\text{Ca}_3(\text{PO}_4)_2$.




$\text{Ca}_3(\text{PO}_4)_2$ is soluble in mineral acids, and in the absence of large excess of the acetate ions, it dissolves in acetic acid too.




2. Carbonate ion solution (use $(\text{NH}_4)_2\text{CO}_3$ or Na_2CO_3 sol.)

The carbonate ion precipitates from slightly basic solution the calcium ion as a white precipitate of calcium carbonate, CaCO_3 .




CaCO_3 is soluble in mineral acids, acetic acid and in a saturated aqueous solution with carbene dioxide.



3. Oxalate ion solution. (use $(\text{NH}_4)_2\text{C}_2\text{O}_4$ sol.)

In ammoniacal solution, the oxalate ion, $\text{C}_2\text{O}_4^{=}$, reacts with the calcium ion to form a white precipitate of calcium oxalate, $\text{CaC}_2\text{O}_4 \cdot \text{H}_2\text{O}$.



4. Flame test:

Calcium salts, especially the chloride, impart when heated in bunsen flame. Impart a brick-red color to the flame.

D . Identification of the Magnesium ion, Mg^{++}

A list of some of the insoluble Mg^{++} compounds is given in table 20


Table 20 Insoluble Magnesium compounds with their colors

Compound	color
* MgC_2O_4	White
MgCO_3	White
MgF_2	White
$\text{Mg}_3(\text{PO}_4)_2$	White
$\text{Mg}(\text{OH})_2$	White
\downarrow $\text{MgNH}_4 \text{ PO}_4$	White


* Decreasing Concentration of Mg^{++} Ion

1. Phosphate ion solution (use $(\text{NH}_4)_2\text{HPO}_4$).

In an ammoniacal solution containing ammonium chloride, magnesium ion reacts with HPO_4^- ion producing a white crystalline precipitates of magnesium ammonium phosphate $\text{MgNH}_4\text{PO}_4 \cdot 6\text{H}_2\text{O}$.




$\text{MgNH}_4\text{PO}_4 \cdot 6\text{H}_2\text{O}$ is soluble in mineral acids, and also in acetic acid too.



2. Aqueous ammonia


Aqueous ammonia precipitates magnesium ions, from solutions containing no ammonia salts, as a white precipitate of magnesium hydroxide, $\text{Mg}(\text{OH})_2$.



The precipitate is readily soluble in acids and in aqueous solutions of ammonia salts.

3. Magneson I (p-nitrobenzeneazoresorcinol) solution.

Sodium or potassium hydroxide precipitates magnesium ion as magnesium hydroxide, $\text{Mg}(\text{OH})_2$



Magneson I, p-nitrobenzeneazoresorcinol, is an organic dye which can be adsorbed by $\text{Mg}(\text{OH})_2$ in a strongly alkaline solution producing a blue precipitate, or rather a blue lake. Other cations of group IV do not interfere in this test.

The principle of separation

It was previously mentioned that the alkaline earth metals and magnesium are included within one group of the periodic table. Thus, they resemble each other closely in their reactions, especially barium, strontium and calcium. This characteristic of close chemical relationship makes it very difficult to obtain a clean separation of these three cations from each other by precipitation methods. This difficulty arises from the fact that any precipitant chosen may bring down all the three cations as precipitates having close values of solubility product constants.

To establish a clean separation between two or more ions using one precipitant, the ratio between the two solubility product of at least one million to one. The ratio of the solubility product constants of any two carbonate or oxalate salts of the last table, is not of the order that permits fairly quantitative separation. Therefore, and because Ca^{++} , Sr^{++} and Ba^{++} ions show characteristic colors of the flame, the flame test which is more specific should be employed.

Slightly soluble carbonate and oxalate salts of group IV cations with their solubility product constants values.


Salt	Ks.p	Salt	Ks.p
MgCO ₃	1.0 x 10 ⁻⁵	MgC ₂ O ₄	8.6 x 10 ⁻⁵
CaCO ₃	8.7 x 10 ⁻⁹	CaC ₂ O ₄	2.3 x 10 ⁻⁹
SrCO ₃	1.6 x 10 ⁻⁹	SrC ₂ O ₄	5.8 x 10 ⁻⁸
BaCO ₃	8.1 x 10 ⁻⁹	BaC ₂ O ₄	1.5 x 10 ⁻⁷

Although several anions can precipitate four members of group IV cations, only the phosphate anion is capable of accomplishing satisfactory precipitation.

For the separation of group IV cations from group V cations, two general methods are known to be used. In the first method, the carbonate anion, CO₃²⁻, is employed as the group precipitant. In this case, the Mg⁺⁺ ion becomes a member of group V cations since it does not precipitate as carbonate. In the second method, the phosphate anion is the precipitate agent for group IV cations.

The precipitation with the phosphate ion takes place in a strongly ammoniacal solution where the calcium, strontium and barium ions are precipitated as the tertiary phosphate while magnesium ion is precipitated as magnesium ammonium phosphate.

Analysis of group IV
Filtrate of group IV and V cations



Procedure NO. 4 Analysis of Group IV cations

<p>Precipitation of the group:</p> <p>The centrifugate (from group III, procedure No.3), may contain any or all the cations of groups IV and V. Adjust the volume to 1 ml discarding any precipitate that may form add 10 drops of 0.5M $(\text{NH}_4)_2\text{HPO}_4$ solution and render strongly basic with concentrated aqueous ammonia. Allow the mixture to stand for a few minutes, test for complete precipitation then centrifuge .</p>		
<p>Precipitate:</p> <p>Contains $\text{Ca}_3(\text{PO}_4)_2$, $\text{Sr}_3(\text{PO}_4)_2$ and $\text{MgNH}_4\text{PO}_4 \cdot 6\text{H}_2\text{O}$ wash with 5 drops of water and discard washings. Add upon the precipitate 5 drops of concentrated CH_3COOH, and stir thoroughly. Any undissolved precipitate will prove by dissolve upon the dilution. Dilute the solution to 2 ml, then add 5 drops of 1M K_2CrO_4 sol. and stir for a minute centrifuge and test for complete pptn.</p>	<p>Centrifugate</p> <p>Contains group V cations.</p>	
<p>Precipitate</p> <p>May contain BaCrO_4 (y). wash with hot water until washings became colorless. Add 5 drops of conc. HCl stir and heat then apply the flame test, a yellowish-green flame confirms Ba^{++} ion.</p>	<p>Centrifugate</p> <p>Contains Ca^{++}, Sr^{++}, Mg^{++}, H_2PO_4^- and CrO_4^{--} ions. Add conc. Aqueous ammonia until the solution is strongly ammoniacal, then centrifuge.</p>	
<p>Precipitate</p> <p>May contains $\text{Ca}_3(\text{PO}_4)_2$, $\text{Sr}_3(\text{PO}_4)_2$ and $\text{MgNH}_4\text{PO}_4 \cdot 6\text{H}_2\text{O}$. dissolve in conc. Acetic acid, as above and dilute to 1 ml, add 5 drops of 0.1M $(\text{NH}_4)_2\text{HPO}_4$ solution and heat to boiling. Test for complete precipitation then cool and centrifuge.</p>	<p>Centrifugate</p> <p>CrO_4 ions Discard</p>	
<p>Precipitate</p> <p>May be SrSO_4 (w). treat with 3 drops of conc. HCl stir and heat then apply the flame test, a crimson-red flame confirms Sr^{++} ion.</p>	<p>Centrifugate</p> <p>May contain Ca^{++} and Mg^{++} ions. add 5 drops of 0.2M $(\text{NH}_4)_2\text{HPO}_4$ solution and heat to boiling. Test for complete precipitation then cool and centrifuge.</p>	
<p>Precipitate</p> <p>May be CaC_2O_4 (w). treat with 2 drops of conc. HCl stir and heat then apply the flame test, a brick-red flame confirms Ca^{++} ion.</p>	<p>Centrifugate</p> <p>May contain Mg^{++} ions. Add 2 drops of magneson I sol. then 3 ml of 3M KOH sol. heat for a few minutes and centrifuge. Wash with 2 ml of KOH solution or water a blue ppt. confirms Mg^{++} ion</p>	

Report sheet

Analysis of group IV cations

Name:.....

Unknown No.:.....

Section:.....

Date:.....

Notes and Equations :

<i>Ions of group IV</i>	<i>Ions found, by the student</i>	<i>Ions present in the unknown solution</i>	<i>grade</i>
Ca ⁺⁺			
Sr ⁺⁺			
Ba ⁺⁺			
Mg ⁺⁺			

Signature of the student

Signature of the grader

Group V Cations

The alkaline metals group
 Na^+ , K^+ and NH_4^+ ions

This group, which may be called soluble group of the cations, is composed of lithium, sodium, potassium and ammonium ions. However, lithium ion will not be included in this manual since it is, known to be a less common compared with sodium and potassium ions.

The alkali metals are the most strongly electropositive elements known, hence they are highly reactive and can not exist in air as free metals due to the rapid reaction with oxygen and moisture, they always tend to attain their stability by losing their outermost shell electron to become ion of +1 oxidation state. They are also characterized by the formation of salts, which, with a few exceptions, are quite soluble in water.

Ammonium ion, NH_4^+ , shows similar properties to those of the alkali metals, such as having a charge of +1, forming salts the majority of which are water soluble and is not precipitated in the previous four groups.

There is no precipitating agent for this group, and identification of each ion is performed in the presence of the other.

A. Identification of the Sodium ion Na^+

A list of some of the slightly soluble sodium salts is given in table 21 below:

Table 21 slightly soluble sodium salts with their colors and molar solubilities

Compound	color	Molar solubility
Na_2SiF_6	White	0.035
$\text{Na}_2\text{H}_2\text{Sb}_2\text{O}_7$	White	0.008

1. Flame test:

Sodium salts, especially the chloride, when heated in bunsen flame. Impart an intense yellow color to the flame.

B. Identification of the Potassium ion K^+

A list of some of the slightly soluble potassium salts is given in table 22 below:

Table 22 slightly soluble potassium salts with their colors and molar solubilities

Compound	color	Molar solubility
KClO_4	White	0.15
K_2PtCl_6	Yellow	0.03
K_2SiF_6	White	0.006
$\text{K}_2\text{CO}(\text{NO}_2)_6$	White	0.001

1. Flame test:

Potassium salts, especially the chloride, when heated in bunsen flame. Impart a pale violate color to the flame. If sodium ions are present, the intense yellow flame of sodium will mask the violate color of potassium flame. Therefore, it will be necessary to view through a filter of two thickness cobalt glass which will absorb the yellow sodium rays leaving only a visible reddish-violate color of the potassium rays.


C. Identification of the Ammonium ion NH_4^+

A list of some of the slightly soluble ammonium salts is given in table 23 below:
Table 23 slightly soluble ammonium salts with their colors and molar solubilities


Compound	color	Molar solubility
$\text{NH}_4\text{HC}_4\text{H}_4\text{O}_6$	White	0.06
$(\text{NH}_4)_2\text{PtCl}_6$	Yellow	0.01
$(\text{NH}_4)_2\text{CO}(\text{NO}_2)_6$	Yellow	0.001

1. Alkali hydroxide solution.

Addition of sodium or potassium hydroxide solution to a solution containing ammonium ions liberates ammonia gas, NH_3 , which can be driven out easily upon gentle heating.




The librated ammonium can be identified by its odor or by its action upon a wet red litmus paper, held at the mouth of the test tube which will turn blue due to the formation of OH^- ions. Another sensitive method is the development of a black color on a filter paper wetted with mercurous nitrate solution, when held on the mouth of the test tube that gives of NH_3 gas.



2. Nessler's reagent (sol. of $\text{K}_2\text{HgI}_4 + \text{NaOH}$)

The addition of Nessler's reagent solution containinh ammonium ions or ammonia produce a brown-red colloidal precipitate or coloration ranging from yellow to brown depending upon the amount of ammonia, or the amount of the ammonium ions present in the solution.



The principle of separation


It is well understood now that the metal ions (cations) of group V are not precipitated with the cations of the earlier groups. In other words, their chlorides, sulfides, hydroxides and phosphate salts are soluble under the conditions of the precipitation by the early group reagents. For this reason, group V is sometimes called the soluble group of the cations.

The test for the NH_4^+ ion must be made on a portion of the original solution, otherwise the test for this cation in group V will always be positive, regardless of its presence in the original solution or not. The test for NH_4^+ ion in group V is meaningless since solution of ammonia or ammonium salts are added several times throughout the group separation. The identification of the ammonium ion by the addition of KOH or NaOH solution is very simple and sensitive, and the evolved NH_3 gas can be recognized by several methods.

Sodium and potassium ions can be identified, in the presence of each other, by means of precipitation reactions. However, the interference of several cations with these tests necessitate the use of the flame test which is a more reliable method for the confirmation of Na^+ and K^+ ions than that of precipitation method.

Analysis of group V cations

Groups V Cations (filtrate from group IV)
 Na^+ , K^+ and NH_4^+ ions



Report sheet

Analysis of group V cations

Name:.....

Unknown No.:.....

Section:.....

Date:.....

Notes and Equations :


<i>Ions of group V</i>	<i>Ions found, by the student</i>	<i>Ions present in the unknown solution</i>	<i>grade</i>
Na ⁺			
K ⁺			
NH ₄ ⁺			

Signature of the student

Signature of the grader

Analytical chemistry

is the science of the characterization and measurement of chemicals and also involve separating ,identifying and determining the relative amounts of the components in a sample of matter. Chemical analysis is divided into two types
Quantitative analysis and *Qualitative analysis* .



Volumetric analysis

Volumetric analysis is a general term for a method in quantitative chemical analysis in which the amount of a substance is determined by the measurement of the volume that the substance occupies. It is commonly used to determine the unknown concentration of a known reactant. Volumetric analysis is often referred to as **titration**.

What is the meaning of Titration?

Titration is a common laboratory method of **quantitative chemical analysis** that is used to determine the unknown **concentration** of a known **reactant**. Because volume measurements play a key role in titration, it is also known as *volumetric analysis*.

A **reagent**, called the *titrant* or *titrator*,^[1] of a known concentration (a **standard solution**) and **volume** is used to react with a solution of the **analyte** or *titrand*,^[2] whose concentration is not known. Using a calibrated **burette** or **chemistry pipetting syringe** to add the titrant.

* A **primary standard solution** is a highly purified compound that serve as a reference material in all volumetric titrimetric methods.

Important requirements for a primary standard are :

- 1-High purify.
- 2-Stability toward air.
- 3-Absence of hydrate water.
- 4-Ready availability at modest cost.
- 5-Reasonable solubility in the titration medium.
- 6-Reasonable large molar mass so that the relative error associated with weighing the standard is minimized.

It is possible to determine the exact amount that has been consumed when the *endpoint* is reached. The **endpoint** is the point at which the titration is complete, as determined by an indicator. This is ideally the same volume as the **equivalence point**—the volume of added titrant at which the number of **moles** of titrant is equal to the number of moles of analyte.

For example:- in the classic strong acid-strong base titration, the endpoint of a titration is the point at which the pH of the reactant is just about equal to 7, and often when the solution takes on a persisting solid color as in the pink of phenolphthalein **indicator**. There are however many different types of titrations.

Many methods can be used to indicate the endpoint of a reaction; titrations often use visual indicators (the reactant mixture changes color). In simple **acid-base titrations** a pH indicator may be used, such as **phenolphthalein**, which becomes pink when a certain pH (about 8.2) is reached or exceeded. Another example is **methyl orange**, which is red in acids and yellow in alkali solutions.

Not every titration requires an indicator. In some cases, either the reactants or the products are strongly colored and can serve as the "indicator". For example, a [redox titration](#) using [potassium permanganate](#) (pink/purple) as the titrant does not require an indicator. When the titrant is reduced, it turns colorless. After the equivalence point, there is excess titrant present. The equivalence point is identified from the first faint persisting pink color (due to an excess of permanganate) in the solution being titrated.

Types of titrations

There are various sorts of titrations whose goals are different to the others. The most common types of titrations in qualitative work are [acid-base titrations](#), [redox titrations](#), complexometric titration and precipitation titration.

1- Acid-base titration

These titrations are based on the neutralization reaction that occurs between an acid and a base, when mixed in solution. The acid is added to a [burette](#) which was rinsed with the same acid prior to this addition to prevent contamination or diluting of the acid being measured. The base is added to a volumetric flask which had been rinsed with distilled water prior to the addition to prevent contamination or dilution of the base/alkali being measured. The solution in the volumetric flask is often a [standard solution](#) (one whose concentration is exactly known). The solution in the burette, however, is the solution whose concentration is to be determined by titration. The indicator used for such an acid-base titration often depends on the nature of the constituents. Common indicators, their colours, and the pH range in which they change colour, are given in the table below, when more precise results are required, or when the titration constituents are a weak acid and a weak base, a [pH meter](#) or a conductance meter are used.

Indicator	Color on Acidic Side	Range of Color Change	Color on Basic Side
Methyl Violet	Yellow	0.0 - 1.6	Violet
Bromophenol Blue	Yellow	3.0 - 4.6	Blue
Methyl Orange	Red	3.1 - 4.4	Yellow
Methyl Red	Red	4.4 - 6.2	Yellow
Litmus	Red	5.0 - 8.0	Blue
Bromothymol Blue	Yellow	6.0 - 7.6	Blue
Phenolphthalein	Colorless	8.3 - 10.0	Pink
Alizarin Yellow	Yellow	10.1 - 12.0	Red

2- Precipitation titration

Precipitation titration :is titration depend upon the combination of ions to form a simple precipitate. Mohr method is a method depend upon formation a colored precipitate for the determination of chloride ion.

A small quantity of potassium chromate (K_2CrO_4) solution is added to serve as indicator. The first excess of titrant results in the formation of a red silver chromate precipitate which signal the end point.

3-Redox titration

These titrations are based on a redox reaction between an oxidizing agent and a reducing agent. The oxidizing agent is added to the burette which was rinsed with the same oxidizing agent. The reducing agent is added to the conical flask, which had been rinsed with distilled water. Like in an acid-base titration, the standard solution is often the one in the conical flask, and the solution whose concentration is to be determined is the one in the burette.

Some redox titrations do not require an indicator, due to the intense colour of some of the constituents. For instance, in a titration where the oxidizing agent **potassium permanganate** (**permanganometry**) is present, a slight faint persisting pink colour signals the endpoint of the titration, and no particular indicator is therefore required.

4- Complexometric titration

These titrations are based on the formation of a **complex** between the analyte and the titrant. The **chelating agent EDTA** is very commonly used to titrate metal ions in solution. These titrations generally require specialized **indicators** that form weaker complexes with the analyte. A common example is **Eriochrome Black T** and murexide for the titration of **calcium** and **magnesium** ions.

Titrimetric analysis generally involve the following steps:

- Sampling
- Titrant preparation
- Standard preparation and conversion to a measurable form
- Titrant standardization by titration of an accurately know quantity of standard
- Sample preparation and conversion to a measurable form
- Sample titration with the titrant solution
- Data analysis

Successful Titrimetric Analysis

A few rules of thumb for designing a successful titration are:

- The titrant should either be a standard or should be standardized.
- The reaction should proceed to a stable and well defined equivalence point.
- The equivalence point must be able to be detected.
- The titrant's and sample's volume or mass must be accurately known.

- The reaction must proceed by a definite chemistry. There should be complicating side reactions.
- The reaction should be nearly complete at the equivalence point. In other words, chemical equilibrium favors products.
- The reaction rate should be fast enough to be practical.
- Substances of known purity for the preparation of standard solution.
- A visual indicator or an instrumental method detecting the completion of the reaction
- Calibrated measuring vessels, including burettes, pipettes, and measuring flasks as follows:-



Burette



Erlenmeyer flask



Volumetric flask

Calculations of volumetric analysis:

Standard solution is one, which contains a known weight of the reagent in a definite volume of the solution.

Molar solution is one, which contains 1 gm molecular weight of the reagent per liter of solution.

$$M = \frac{\text{Weight}}{\text{M.Wt}} * \frac{1000}{\text{Volume(mL)}}$$

Normal solution is one that contains 1gm equivalent weight per liter of solution.

$$N = \frac{\text{Weight}}{\text{eq.wt}} * \frac{1000}{\text{Volume(mL)}}$$

Part per million(ppm): Milligrams of solute per liter of solution.


$$\text{ppm} = \frac{\text{weight solute(mg)}}{\text{Volume solution(L)}}$$

$$\text{ppm} = \frac{\text{wt(g)}}{\text{V(mL)}} \times 10^6$$

$$\text{ppm} = M \times \text{M.wt} \times 1000$$

.....

For titrimetric reaction:



At equivalent point:

no.mmol of titrant(A)= no.mmol of titrand(B)

$$N_A \times V_A = N_B \times V_B$$

$$\text{or } M_A \times V_A \times R = M_B \times V_B \quad (R=b/a)$$

Equivalent weights

(1)Equivalent weight in neutralization reactions.

The equivalent weight of acid is that weight of it which contains one-gram atom of replaceable hydrogen.

Ex: equivalent weight of $\text{H}_2\text{SO}_4 = \text{M.Wt H}_2\text{SO}_4/2$
equivalent weight of $\text{H}_3\text{PO}_4 = \text{M.Wt H}_3\text{PO}_4/3$

$\text{eq.wt acid} = \frac{\text{M.Wt acid}}{\text{No. of active H}^+}$

The equivalent weight of Base is that weight of it which contains one replaceable hydroxyl group.

Ex: equivalent weight of $\text{NaOH} = \text{M.Wt NaOH}/1$


$\text{eq.wt base} = \frac{\text{M.Wt base}}{\text{No. of active OH}^-}$
--

(2)Equivalent weight in Oxidation -reduction reactions.

The equivalent weight of an oxidant or a reductant is the number of electrons which 1mol of the substance gains or loses in the reaction.

$\text{eq.wt} = \frac{\text{Formula weight (gm/mole)}}{\text{No. of lost or gained electrons}}$

Ex:




$$\text{eq.wt MnO}_4^- = \frac{\text{F.wt MnO}_4^-}{5} \qquad \text{eq.wt C}_2\text{O}_4^{2-} = \frac{\text{F.wt C}_2\text{O}_4^{2-}}{2}$$

(3) Equivalent weight of complex formation and precipitation reactions.

Here the equivalent weight is the weight of the substance which *contains* or *reacts* with 1g. atm of a univalent cation M^+ .

Ex: When silver nitrate reacts with sodium chloride, to form silver chloride, the equivalent weight of AgNO_3 is:



.....

Experiment (1)

A: Preparation of approximately (0.1N) HCl

1- Calculate the normality of the concentrated HCl :

$$N = \frac{\text{Sp.g} * \text{Percentage} * 1000}{\text{equivalent weight}}$$

$$N = \frac{1.19 * (37/100) * 1000}{36.5} = 12.0630$$

2-To prepare (500mL) of 0.1N HCl :

$$N_1 \times V_1 = N_2 \times V_2$$

conc.HCl dil.HCl

$$12.0630 \times V_1 = 500 \times 0.1$$
$$V_1 = 4.1449 \text{ mL}$$

Transfer distilled water to a glass-Stoppard bottle and add about 4.14 ml of conc. HCl Mix the contents of the bottle thoroughly. Dilute to the mark with distilled water and make a label.

B: Preparation of (0.1N) sodium carbonate:

1-Dry a quantity of primary-standard sodium carbonate for 2 hours at 110 °C, cool in a desiccator .Weigh 1.325gm of the salt into a small beaker, and dissolve in 25 to 50 mL of distilled water.

2-Transfer this solution to a 250mL volumetric flask. Wash the beaker and transfer this washing water to the volumetric flask. Dilute to the mark with distilled water and mix thoroughly. Make a label.

$$N = \frac{\text{Weight}}{\text{equivalent weight}} * \frac{1000}{\text{Volume(mL)}}$$

Experiment (2)

Standardization of HCl solution with standard solution of Na₂CO₃

- 1-Clean the burette and rinse with HCl solution.
- 2-Fill the burette with HCl.
- 3-Pipet 10 ml of standard solution (Na₂CO₃) into a 250mL conical flask. Add 3 drops of Methyl orange indicator.
- 4- Titrate by adding HCl drop wise until the solution just begins to change from yellow to red.
- 5-Repeat the titration a few times until you get approximate results. Take the average of the results and subtract 0.05mL. (This result represents the volume of extra drop which change the color of indicator).
- 6- Calculate the normality of HCl :
$$N_{\text{acid}} \times V_{\text{acid}} = N_{\text{base}} \times V_{\text{base}}$$
- 7- Make label on your bottle containing your name, date of preparation and concentration of acid after standardization .
- 8- The equation of reaction :



Experiment (3)

Analysis of sodium carbonate Na₂CO₃

- 1-Clean the burette and rinse with standardized HCl solution and then fill it with the acid.
- 2-Pipet 10 ml unknown solution (Na₂CO₃) into a 250mL conical flask . Add 2 drops of phenolphthalein indicator the solution will be pink.
- 3- Titrate by adding HCl drop wise until the solution just begins to change its color from pink to colorless this data will be (V1).
- 4- Add 1-2 drops of Methyl orange indicator to the above solution which became yellow then complete the titration until the color of the solution became pale orange (onion), this data will be (V2).
- 5-Repeat the titration a twice time until you get approximate results. Take the average of the results and subtract 0.05mL. (This result represents the volume of extra drop which change the color of indicator).
- 6- Make a table as bellow:-

	First titration	Second titration	Third titration	Average of titrations
Titration with ph.ph	V1	V1	V1	V1 (av.)

Titration with M.O	V2	V2	V2	V2 (av.)
-----------------------	----	----	----	----------

$$V1(av.) = \frac{V1+V1+V1}{3}$$


$$V1(av.) = 1/2 CO_3^{=}$$

$$V2(av.) = 1/2 CO_3^{=}$$

$$V1+V2 = V_{tot.} \text{ of } Na_2CO_3$$

$$N_{acid} \times V_{(tot.) acid (from burette)} = N_{base} \times V_{base}$$

7- The equation of reaction :



Experiment (4)

Analysis of a mixture (sodium hydroxide + sodium carbonate)

1-Transfer a (10ml) of a mixture solution to 250 ml conical flask and add 1drop of phenolphthalein indicator.

2-Clean the burette and rinse with hydrochloric acid.

3-Fill the burette with HCl.

4-Titrate with standard hydrochloric acid solution until the pink color disappears (the solution will be colorless).

* At this stage: all the hydroxide and half the carbonate have been neutralized. Let us assume that the volume of acid be "X" ml.

5-Now add 3drops of Methyl Orange indicator into the solution above and continue the titration until the solution just begins to change from yellow to red.

* At this stage :another half of carbonate has been neutralized. Let the volume of acid be "Y" ml.

6-Repeat the titration a few times until you get approximate results.

Calculations:

$$(X-Y)=Z \rightarrow \text{Volume of HCl which equivalent to } OH^- \dots\dots\dots(1)$$

$$2 \times Y \rightarrow \text{Volume of HCl which equivalent to } CO_3^{=}\dots\dots\dots(2)$$

$$N_{acid} \times V_{acid} = N_{base} \times V_{base}$$

$$N_{acid} \times Z = N_{OH^-} \times 10 \dots\dots\dots (1)$$

$$N_{acid} \times 2Y = N_{CO_3^{=}} \times 10 \dots\dots\dots (2)$$

$$\text{Concentration of } [OH^-]_{(ppm)} = N_{OH^-} \times \text{eq.wt} \times 1000$$

$$\text{Concentration of } [\text{CO}_3^{2-}]_{(\text{ppm})} = N \text{ CO}_3^{2-} \times \text{eq.wt} \times 1000$$

.....

Experiment (5)

Analysis of a mixture (sodium bicarbonate + sodium carbonate)

1-Repeat the same steps of the experiment no.3.

2-Volume of acid at the end point of Ph.Ph indicator=X.

3-Volume of acid at the end point of M.O indicator=Y.

X=Volume of HCl \equiv **1/2 CO₃²⁻**

Y=Volume of HCl \equiv **1/2 CO₃²⁻ + HCO₃⁻**

(Y-X)=Z \rightarrow Volume of HCl which equivalent to HCO₃⁻(1)

2×X \rightarrow Volume of HCl which equivalent to CO₃²⁻..... (2)

$$N_{\text{acid}} \times V_{\text{acid}} = N_{\text{base}} \times V_{\text{base}}$$

$$N_{\text{acid}} \times Z = N \text{ HCO}_3^- \times 10 \quad \text{..... (1)}$$

$$N_{\text{acid}} \times 2X = N \text{ CO}_3^{2-} \times 10 \quad \text{..... (2)}$$

Concentration of [HCO₃⁻]_(ppm) = N HCO₃⁻ × eq.wt × 1000

Concentration of [CO₃²⁻]_(ppm) = N CO₃²⁻ × eq.wt × 1000

.....

Experiment (6)

Determination of chloride ion by Mohr method

Precipitation titration :is titration depend upon the combination of ions to form a simple precipitate. Mohr method is a method depend upon formation a colored precipitate for the determination of chloride ion.

Chloride ion ,is reacted with silver nitrate solution to form AgCl precipitate.



$$K_{\text{sp AgCl}} = 1.2 \times 10^{-10}$$

A small quantity of potassium chromate (K₂CrO₄) solution is added to serve as indicator. The first excess of titrant results in the formation of a red silver chromate precipitate which signal the end point.



$$K_{\text{sp Ag}_2\text{CrO}_4} = 1.7 \times 10^{-12}$$

Procedure:

1-Clean the burette and fill it with silver nitrate(0.1N).

2-Pipet 10mL of chloride ion solution into 250mL conical flask, add 5drops of potassium chromate.

3-Titrate chloride solution against silver nitrate until arrive to equivalent point(the point in which the number of moles of AgNO₃ equal to the number of moles of chloride ion)[Notice a white precipitate in the yellow solution].

After this point the excess of AgNO₃ will react with potassium chromate leading to formation of red precipitate Ag₂CrO₄ (the end point).The difference between equivalent and end point is the volume of AgNO₃ reacted with the indicator.

4- Repeat the titration and calculate the average and subtract the indicator blank(equal to 0.1mL).Calculate the normality.

$$N_{Cl^-} \times V_{Cl^-} = N_{Ag^+} \times V_{Ag^+}$$

$$\text{Concentration of } [Cl^-]_{(ppm)} = N_{Cl^-} \times \text{eq.wt} \times 1000$$

Experiment (7)

Oxidation-reduction reactions

Oxidation-reduction reactions: are all the reactions involving change in oxidation number or transfer of electrons among the reacting substances.

* **Reduction agent:** is one that loses electrons and is oxidized to a higher valiancy conditions.

Ex: Ferrous, stannous compounds, sodium thiosulphate and oxalate ion.



* **Oxidation agent:** is one that gains electrons and is reduced to a lower valiancy conditions.

Ex: Potassium dichromate, potassium permanganate.

The permanganate ion in acid solution is a strong oxidizing agent, because the standard potential in acid solution E° has been calculated to be 1.51V.



Preparation of solutions:

A: Preparation of 0.1N potassium permanganate.

* Weigh out about 3.2gm of $KMnO_4$ on a watch glass ,and transfer it to a 1500mL beaker and add 1L distilled water ,cover the beaker with a clock glass then heat the solution to boiling, boil gently for 1hr and allow the solution to cool to the laboratory temperature , filter the solution through a funnel containing a plug of purified glass wool, collect the filtrate in a vessel. The solution should be stored in a clean ,glass stoppered bottle of a dark-brown-colored glass and kept in the dark.

B: Preparation of 0.1N sodium oxalate ($Na_2C_2O_4$).

* Weigh out about 1.675gm of $Na_2C_2O_4$ and dissolved it in distilled water in a 250mL volumetric flask.

$$N = \frac{Wt}{\text{eq.wt}} * \frac{1000}{V(\text{mL})}$$

$$0.1 = \frac{Wt}{67} * \frac{1000}{250} \implies Wt = 1.675 \text{ gm}$$

Standardization of permanganate solution with oxalate ion:

1-Pipet 10mL of 0.1N of sodium oxalate ($Na_2C_2O_4$) in a 250mL conical flask .




2-Add by a graduated cylinder 5mL of H_2SO_4 [1:8].

3-Add some drops of potassium permanganate and heat to 55-60°C and complete the titration by adding permanganate solution until a faint pink color persists for 30sec.

4-Calculate the normality of permanganate ion.

$$N_1 \times V_1 = N_2 \times V_2$$

$$\text{KMnO}_4 \quad \text{KMnO}_4 \quad \text{Na}_2\text{C}_2\text{O}_4 \quad \text{Na}_2\text{C}_2\text{O}_4$$



Determination the concentration of ferrous ion:



1-Pipet 10mL of ferrous ion solution in a 250mL conical flask.

2-Add 5mL of sulphuric acid H_2SO_4 [1:8].

3-Titrate the solution with standard (0.1N) KMnO_4 until the pink color persists.

4-Repeat the titration and then calculate the normality of ferrous ion.

$$V_{\text{KMnO}_4} \times N_{\text{KMnO}_4} = \frac{\text{Wt Fe}^{2+}}{\frac{55.85}{1000}}$$



.....

Experiment (8)

Complexometric titration

This kind of titration is called complexometric because it involved formation of complex especially by using the compound EDTA (Ethylene Diamine Tetra Acetic acid) which forms a stable complex with metal ions especially with calcium and magnesium ions. In addition the users indicator with EDTA have also capable to form a complex with the same ions (Mg^{2+} and Ca^{2+}) but less stable than the complex formed between the metal ions and EDTA.

The total hardness of water is generally due to dissolved calcium and magnesium salts. Water hardness was defined in terms of the capacity of cations in the water to replace the sodium and potassium ions in soaps and form sparingly soluble products.

There are two types of water hardness, temporary and permanent.

Temporary hardness is due to the bicarbonate ion, HCO_3^- , (sodium and magnesium bicarbonate) being present in the water. This type of hardness can be removed by boiling the water to expel the CO_2 .

Permanent hardness is due to the presence of the ions Ca^{2+} , Mg^{2+} and SO_4^{2-} . This type of hardness cannot be eliminated by boiling.

Hardness is now expressed in terms of the concentration of calcium carbonate that is equivalent to the total concentration of all the multivalent cations in the sample. The determination of hardness is a useful analytical test that provides a measure of water quality for household and industries uses.

Determination of total hardness (permanent and temporary) of water

A: Determination of total hardness of water.

- 1-Prepare 250mL of 0.02N of EDTA solution.
- 2-Clean the burette and rinse with EDTA solution.
- 3-Pipet 10mL of hardness water and add 2mL of ammonia buffer solution (PH=10)[prepared from ammonium chloride +conc. ammonia].
- 4-Add 3-4 drops of Eriochrome black-T indicator and titrate with standard EDTA solution until the color changes from wine red to pure blue color.
- 5-Repeat the titration a few times and calculate the total hardness ($V_{\square A}$) of the sample of water ($\text{Mg}^{2+} + \text{Ca}^{2+}$).

Calculations

$V_{\text{EDTA}} \equiv \text{Volume of EDTA which equivalent to } (\text{Mg}^{2+} + \text{Ca}^{2+}) \Rightarrow V_{\square A}$

$$\begin{array}{ccccccc} N_1 & \times & V_1 & = & N_2 & \times & V_2 \\ \downarrow & & \downarrow & & \downarrow & & \downarrow \end{array}$$

EDTA EDTA Sample(?) Sample
(0.02N)

Total hardness (ppm) = $N \times \text{eq. wt} \times 1000$
(Concentration of $\text{Mg}^{2+} + \text{Ca}^{2+}$) CaCO_3

Or **Total hardness (ppm) = $(V_{\text{EDTA}}/V_{\text{solution}}) \times 2000$**
(Concentration of $\text{Mg}^{2+} + \text{Ca}^{2+}$)

B: Determination of the permanent hardness

- 1-Pipet 10mL of the hardness water and add 2mL of sodium hydroxide solution 4N (PH=12) which is sufficiently high to precipitate the magnesium quantitatively as the hydroxide.
- 2-Add some crystals of Murexide indicator.
- 3-Stirr the solution and titrate it with standard EDTA solution until the color changes from pink to a purple.
- 4-Repeat the titration a few times and calculate the ($V_{\square B}$), which is the volume of EDTA equivalent to Ca^{2+} only.

Calculations

$V_{\text{EDTA}} = \text{Volume of EDTA which equivalent to } \text{Ca}^{2+} \text{ only} \Rightarrow V_{\square B}$

$$N_1 \times V_1 = N_2 \times V_2$$

↓ ↓ ↓ ↓

EDTA EDTA Sample(?) Sample
(0.02N)

Permanent hardness (ppm)=N × eq.wt × 1000
(Concentration of Ca²⁺) CaCO₃

Or Permanent hardness (ppm)=(V_{EDTA}/V_{solution})× 2000
(Concentration of Ca²⁺)

C: Determination of temporary hardness

Calculate the temporary hardness of water by subtracting the permanent hardness from the total hardness.

Temporary hardness = Total hardness – Permanent hardness = V_A - V_B
(Mg²⁺) (Mg²⁺ + Ca²⁺) (Ca²⁺)

N₁ × (V_A - V_B) = N₂ × V₂
 ↓ ↓ ↓ ↓
EDTA EDTA Sample(?) Sample
(0.02N)

Temporary hardness (ppm)=N × eq.wt × 1000
(Concentration of Mg²⁺) CaCO₃

Or Temporary hardness (ppm)=(V_{EDTA}/V_{solution})× 2000
(Concentration of Mg²⁺)

University of Baghdad
College of Science
Department of Chemistry

ANALYTICAL CHEMISTRY LABORATORY

Analytical Chemistry Sheet
For The First Class

Quantitative Analysis

Volumetric Analysis

The second season
2020-2021

